



Innovation • Value • Discovery ● ●



REF K-3245



Innovation • Value • Discovery ● ● ●

ExiPrep™ Bacteria Genomic DNA Kit

ExiPrep™ Bacteria Genomic DNA Kit

Kit for the extraction of bacteria genomic DNA in various samples using ExiPrep™ instrument

User's Manual



Version No.: 3.0 (2010-12-24)

BIONEER CORPORATION



49-3 Munpyeong-Dong, Daedeok-Gu, Daejeon, 306-220, Korea

Tel: +82-42-930-8777

Fax: +82-42-930-8600

Email: order@bioneer.co.kr

www.bioneer.com

Safety Warnings and Precautions

Please inquire BIONEER's Customer Service Center to obtain a copy of the Material Safety Data Sheet (MSDS) for this product.

Before, during and after use of this kit as described in this User Manual, all potentially hazardous materials (i.e. materials that may have come in contact with clinical samples) including tubes, tips and materials should be processed and disposed of according to applicable and appropriate regulations of the municipality/government in which this product is being used.

Please read the User Manual before using this Kit. Please check the integrity of all tubes, tips and other materials supplied with this kit prior to use. Adhere to general clinical laboratory safety procedures during the experiment.

Warranty and Liability

All BIONEER products are manufactured and tested under strict quality control protocols certified under **ISO 9001** and **ISO 13485** standards. BIONEER guarantees the quality of all directly manufactured products during the warranty period of one (1) year from date of purchase. If any issues are discovered relating to compromise in product quality, immediately contact BIONEER's Customer Service Center (order@bioneer.com).

Trademarks

ExiPrep™ is the trademark of Bioneer Corporation. KOREA

Copyright © 2011 by Bioneer Corporation. All rights reserved.

Table of Contents

1.	KIT COMPONENTS	1
2.	STORAGE	2
3.	INTRODUCTION	3
4.	STARTING SAMPLE&TYPICAL YIELD.....	4
5.	GENOMIC DNA EXTRACTION	4
5.1	Additionally required materials & devices.....	4
5.2	Extraction from cultured bacteria.....	5
5.3	Extraction from urine.....	10
5.4	Extraction from sputum	11
5.5	Extraction from swab.....	13
6.	TROUBLESHOOTING	14
7.	PROTOCOL NUMBER LIST.....	16
8.	ORDERING INFORMATION	17
9.	REFERENCES	17
10.	EXPLANATION OF SYMBOLS	18

1. Kit components

	ExiPrep™ Bacteria Genomic DNA Kit (K-3245)
Buffer Cartridge ①	6 ea
Buffer Cartridge ②	6 ea
Resuspension Buffer	1 ea
Disposable Tip	3 pack
Elution Tube (8-strip)	1 pack
User's Manual	1 ea
One Page Protocol	1 ea

2. Storage

ExiPrep™ Kits are designed to contain all necessary buffers (including binding/washing/elution buffers) and components within the Buffer Cartridges. All Buffer Cartridges are sealed with protective film to prevent leakage, evaporation and contamination. The Buffer Cartridges can be stored at room temperature (15°C-25°C) under dry conditions for up to 2 years without opening.

ExiPrep™ Kits may also provide enzymes in lyophilized form for increased convenience.


Lyophilized enzymes may be contained within the Buffer Cartridges and/or 2.0mL screw cap tubes. The enzymes in lyophilized form can be stored at room temperature (15°C-25°C) up to 2 years without any detectable loss in activity. Reconstituted enzymes must be stored at -20°C

The disposable tips and elution tubes included within the Kit are DNase and RNase free. Please take caution in storing the consumables and avoid nuclease contamination.

3. Introduction

ExiPrep™ Bacteria Genomic DNA Kit is designed for the extraction of genomic DNA from gram negative and positive bacteria and fungi. Gram positive bacteria and fungi must undergo an enzymatic digestion step for optimal extraction. The digested samples then must be resuspended with the provided Resuspension Buffer before proceeding with genomic DNA extraction.

Samples are mixed with Lysis Buffer contained within Buffer Cartridge ① and undergo incubation so the nucleic acids may be released into the solution. After the incubation step, silica magnetic beads will be introduced to bind the DNA. The magnetic beads can be separated from the solution (containing cellular debris) through magnetization. The beads are then washed to flush any impurities that may exist. Finally, a small volume of elution buffer washed over the magnetic beads will release the DNA. The eluted DNA is in a highly pure form and is ready to be used.



[Description of DNA extraction using silica magnetic bead]

4. Starting Sample & Typical yield

Typical amounts of starting material quantity, elution volume and the typical yields are described below. Results of your application may vary.

Sample type	Starting Sample	Elution Volume	Typical Yield
Gram (-) bacteria	~1x10 ⁸ cells	50~200ul	5-15ug
Gram (+) bacteria	~1x10 ⁸ cells	50~200ul	5-15ug
Yeast (<i>S. pombe</i>)	~1x10 ⁸ cells	50~200ul	5-15ug
Fungi	~1x10 ⁸ cells	50~200ul	5-15ug



5. Genomic DNA Extraction





5.1 Additionally required materials & devices

- Disposable gloves
- Pipettes
- Sterilized pipette tips with filters
- 1.5ml microcentrifuge tubes or 15ml conical tubes
- 1x PBS buffer, 1x TE buffer, 10N NaOH
- Lysozyme or lyticase
- Table top centrifuge or swing rotor centrifuge
- ExiPrep™ 16 Plus (Cat. No. : A-5030, Bioneer Co.,KOREA),
ExiPrep™ 16 Pro (Cat. No. : A-5040, Bioneer Co.,KOREA)





5.2 Extraction from cultured bacteria

This protocol is designed for extraction of genomic DNA from cultured bacteria.

	<ol style="list-style-type: none"> a. For Gram (-) bacteria Resuspend the cell pellet (up to 1×10^8 cells) in 200ul of Resuspension Buffer. b. For Gram (+) bacteria <ul style="list-style-type: none"> - Resuspend the cell pellet (up to 1×10^8 cells) in 200ul of 1x TE buffer - Add 20ul of lysozyme (50mg/ml) and incubate the tube at 37°C for at least 1hr. - Centrifuge the tube at 5000rpm for 3min and remove the supernatant. - Re-suspend the pellet with 200ul of Resuspension Buffer.
	<ol style="list-style-type: none"> Place the elution tube rack, tip rack, and Buffer Cartridge ① on the setup tray.







	<ol style="list-style-type: none"> Ensure that all tips and tubes are inserted in the desired positions.
	<ol style="list-style-type: none"> Punch holes through the Buffer Cartridge sealing film with the 6-hole punch tool in the locations corresponding to the tubes and tips. <ul style="list-style-type: none"> ※ Before punching the holes, make sure to shake Buffer Cartridge ① to suspend the beads, then gently swing outward to settle solutions to the bottom of the wells.
	<ol style="list-style-type: none"> Load 200 µL of sample into the sample loading wells.
	<ol style="list-style-type: none"> Turn on the ExiPrep™ instrument Press the 'Store' button for cooling. (ExiPrep™ 16 Pro only) Press the 'Start' button to access the PREP SETUP menu.

ExiPrep™ Bacteria Genomic DNA Kit

	<p>9. Enter a protocol number from the protocol list (Chapter 7) that best suits your sample source and target nucleic acid type.</p> <p>10. Press the 'Enter' button to proceed to the next step.</p>
	<p>11. Select the desired elution volume from the screen.</p> <p>※ Recommended elution volume: 100ul</p> <p>12. Press the 'ok' button to move to the next step.</p>
	<p>13. Open the instrument door and pull out the Base Plate.</p> <p>14. Place all racks and Buffer Cartridges in their correct positions on the Base Plate according to the CHECK LIST (steps 15~19).</p>
	<p>15. Place Buffer Cartridge ② into its labeled position on the Base Plate.</p>




7

ExiPrep™ Bacteria Genomic DNA Kit

	<p>16. Place Buffer Cartridge ① into its labeled position on the Base Plate.</p>
	<p>17. Place the Elution Tube Rack into its labeled position on the Base Plate.</p> <p>a. For ExiPrep™ 16 Plus Place the Elution Tubes in Lanes A and B of the Elution Tube Rack.</p>
	<p>b. For ExiPrep™ 16 Pro Place the Elution Tubes into the 2nd and 3rd rows of the Elution Tube Rack.</p>
	<p>18. Place the Disposable Filter Tip Rack into its labeled position on the Base Plate.</p>
	<p>19. Place the Waste Tray into its proper position (between Buffer Cartridges ① and ②).</p>
	<p>20. Push the Base Plate completely back into the instrument and close the door.</p>

8



ExiPrep™ Bacteria Genomic DNA Kit

	<p>21. Verify the protocol name on the screen. The first two letters represent the target nucleic acid type, and the next two letters represent the sample source.</p> <p>22. Press the 'Run' button to start the extraction run.</p>
	<p>23. When the run process is successfully completed, three options will display.</p>
	<p>24. Open the door and pull out the Base Plate.</p> <p>25. Take the Elution Tubes from Base Plate.</p> <p>26. Remove the Buffer Cartridges and all racks from Base Plate and close the door.</p> <p>27. Quantify the purified genomic DNA (if desired).</p>

ExiPrep™ Bacteria Genomic DNA Kit



5.3 Extraction from urine



This protocol is designed for the extraction of genomic DNA from urine samples.

	<ol style="list-style-type: none"> 1. Transfer an adequate volume of urine into a 15ml conical tube. 2. Centrifuge at 3,000 rpm for 20 min. ※ Please repeat centrifugation if pelleting is insufficient. 3. Add 200µl of Resuspension Buffer and mix well by vortexing for 15 sec. ※If the sample has excessive sedimentation, the volume of Resuspension Buffer may be increased up to 500µl.
	<ol style="list-style-type: none"> 4. Go to step 2. of '5.2 Extraction from cultured bacteria'. <p>※ Tip: These steps may be prepared during the sample centrifugation step.</p>

5.4 Extraction from sputum




This protocol is designed for the extraction of genomic DNA from sputum samples.

	<ol style="list-style-type: none"> 1. Add 1/10 volume of 10N NaOH into the container holding the sputum sample and mix by vortexing for 3 min. Incubate the sample tubes for 15 min at room temperature. ※ If sample has high viscosity, add more 10N NaOH (up to 500μl) or incubate the sample tube for up to 30 min.
	<ol style="list-style-type: none"> 2. Transfer 1.5ml of sample into a 1.5ml microcentrifuge tube. 3. Centrifuge the tube at 9,000 rpm for 5 min and completely remove the supernatant. ※ If your centrifugation is incomplete, please repeat centrifugation.

	<ol style="list-style-type: none"> 4. Re-suspend the pellet with 1ml of 1X PBS buffer. 5. Centrifuge the tube at 9,000 rpm for 5 min and completely remove the supernatant. ※ 6. Repeat step 4.~5. over two times. 7. Resuspend the pellet with 200μl of Resuspension Buffer.
	<ol style="list-style-type: none"> 8. Go to step 2. of '5.2 Extraction from cultured bacteria'.

5.5 Extraction from swab

This protocol is designed for the extraction of genomic DNA from swab samples.

	<ol style="list-style-type: none"> 1. Transfer the swab into a 15ml conical tube or a specialized sample tube. 2. Add 1ml of 1x PBS buffer or normal saline and mix well by vortexing. ※ If the specialized sample tube contains media, do not add the 1ml liquid and proceed to the next step.
	<ol style="list-style-type: none"> 3. Transfer 500µl of sample into a 1.5ml microcentrifuge tube. 4. Centrifuge the tube at 13,000 rpm for 2 min and remove the supernatant. 5. Resuspend the pellet with 200µl of Resuspension Buffer.
	<ol style="list-style-type: none"> 6. Go to step 2. of '5.2 Extraction from cultured bacteria'.

6. Troubleshooting

A. Low yield of Genomic DNA

- 1) Did you add sufficient amount of sample? The yield is dependent on the sample type and amount. Overloading sample may sometimes decrease the yield.
- 2) Did you lyse the samples completely? Incomplete lysis decreases yield and purity.
- 3) Is there salt precipitate in any of the buffers? Warm the bottles at 60°C to re-dissolve.
- 4) Did you shake Buffer cartridge ① before use? Incomplete resuspension of the magnetic beads may decrease the yield and purity.

B. Co-eluted magnetic particles

The magnetic particles may co-elute with the extracted DNA after genomic DNA extraction. Co-eluted magnetic particles cannot bind nucleic acid in Elution Buffer and it will not decrease the yield and purity.

Co-eluted magnetic particles can easily be separated by simple centrifugation.

C. Empty elution tubes

1) Too much starting sample

Please add the volume of sample suggested in this manual. Excessive sample volume may cause clogging by magnetic beads and interfere with normal extraction.

2) Insufficient elimination of the supernatant

When processing samples to prepare them for extraction, you must ensure complete removal of sample supernatant (ex. Lipid layer, medium). Insufficient removal of impurities may cause clogging by magnetic beads and interfere with normal extraction.

D. Expansion of sealing film

Please store the Buffer Cartridges in recommended storage conditions (Chapter 2. Storage). High storage temperatures will increase the vapor pressure of certain buffers and force the film to bulge.

7. Protocol number list

No.	Target	Sample source
1 01	Genomic DNA	Whole blood
1 02	Genomic DNA	Animal tissue
1 03	Genomic DNA	FFPE tissue
1 04	Genomic DNA	Plant tissue
1 05	Genomic DNA	Plant seed
1 06	Genomic DNA	Rice
1 07	Genomic DNA	Cultured cell
1 08	Genomic DNA	Gram (+) bacteria
1 09	Genomic DNA	Gram (-) bacteria
1 10	Genomic DNA	Yeast
1 11	Genomic DNA	Fungi
1 14	Genomic DNA	Buffy coat
1 15	Genomic DNA	Sputum
1 16	Genomic DNA	BAL
1 17	Genomic DNA	Saliva
1 18	Genomic DNA	Swab
1 19	Genomic DNA	Urine
1 20	Genomic DNA	Stool
1 23	Genomic DNA	CSF
1 24	Genomic DNA	EPS
1 25	Genomic DNA	Respiratory sample

8. Ordering information

Product	Size	Cat. No.
ExiPrep™ 16Plus	1 ea	A-5030
ExiPrep™ 16Pro	1 ea	A-5040
ExiPrep™ Blood Genomic DNA Kit	96 prep.	K-3215
ExiPrep™ Tissue Genomic DNA Kit	96 prep.	K-3225
ExiPrep™ Bacteria Genomic DNA Kit	96 prep.	K-3245
ExiPrep™ Plant Genomic DNA Kit	96 prep.	K-3255
ExiPrep™ Viral DNA Kit	96 prep.	K-3515
ExiPrep™ Viral RNA Kit	96 prep.	K-3525
ExiPrep™ Beef Genomic DNA Kit	96 prep.	K-3200-CB
ExiPrep™ Rice Genomic DNA Kit	96 prep.	K-3200-CR

9. References

Obata, K., Segawa, O., Yakabe, M., Ishida, Y., Kuroita, T., Ikeda, K., Kawakami, B., Kawamura, Y., Yohda, M., Matsunaga, T. and Tajima, H. (2001). Development of a novel method for operating magnetic particles magtration technology, and its use for automating nucleic acid purification. *J. Biosci. Bioeng.*, 91(5), 500–503.

Honore-Bouakline, S., Vincensini, J. P., Giacuzzo, V., Lagrange, P. H., Herrmann, J. L. (2003). Rapid Diagnosis of Extrapulmonary Tuberculosis by PCR: Impact of Sample Preparation and DNA Extraction. *J. Clin. Microbiol.* 41: 2323-2329

Boom, R., Sol, C., Beld, M., Weel, J., Goudsmit, J. and Wertheim-Van Dillen, P. (1999) Improved silica-guanidiniumthiocyanate DNA isolation procedure based on selective binding of bovine alphacasein to silica particles. *Journal of Clinical Microbiology* 37, 615–619

10. Explanation of symbols



Catalog Number



Batch code



Contains Sufficient for test



Expiration Date



Storage Conditions (Temp.)



Manufacturer



Caution, Consult accompanying documents

● Bioneer World Wide

Headquarters

Address 49-3, Munpyeong-dong, Daedeock-gu, daejeon 306-220, Korea
Tel +82-42-930-8777 (Korea only: 1588-9788)
Fax +82-42-930-8600
E-mail sales@bioneer.co.kr
Web <http://www.bioneer.co.kr>

Seoul Office

Address 2F Samhwawa Sansu Building, 6, Yangjae-Dong, Seocho-Gu, Seoul 137-130, Korea
Tel +82-2-598-1094 (korea only: 1588-9788)
Fax +82-8-598-1096
E-mail sales@bioneer.co.kr
Web <http://www.bioneer.co.kr>

Bioneer, Inc.

Address 1000 Atlantic Avenue Alameda, CA 94501 USA
Tel +1-877-264-4300(Toll Free)
Fax +1-510-865-0350
E-mail order.usa@bioneer.com
Web www.bioneer.com

Bioneer Trade (Shanghai) Co, Ltd

Address 403Room, Building 88, no.887, Zuchongzhi Road, Zhangjiang High Tech Park,
PuDong New District, Shanghai 201203, China
Tel +86-21-5080-0969
Fax +86-21-5080-1620
E-mail asan@bioneer.com
Web www.bioneer.com