

AccuPower® GreenStar™ qPCR PreMix

I. Introduction

AccuPower® GreenStar™ qPCR PreMix is a ready-to-use reagent containing all components for real-time PCR reaction, except for target-specific primers. Just addition of specific primers and target gene into tube provide reproducible results with high sensitivity and specificity. Because all components for PCR reaction are lyophilized in real-time PCR plates or tubes, the stability of product can be extremely extended up to 2 years at -20 °C storage, compared to that of other commercially available product.

This product can be used in real-time PCR experiments for the amplification and detection of genomic DNA and cDNA targets, differential gene expression profiling, and Microbial & Viral pathogen detection. This product provides the reproducible results with the superior specificity, high sensitivity, wide dynamic range and accurate quantification.

II. Principle

PCR products are detected with SYBR Green I dye in real-time monitoring.


1) PCR (Polymerase Chain Reaction)

PCR is a biochemistry and molecular biology technique for amplification of target DNA across several orders of magnitudes, generating millions or more copies of target DNA pieces.

There are three major steps at different temperatures in a PCR, which are repeated for 30 or 45 cycles. Double-stranded target DNA is heat-denatured (denaturation step), the two primers complementary to the target segment are annealed at low temperature (annealing step), and the annealed primers are then extended at an intermediate temperature (extension step) with a DNA polymerase. As the target copy number doubles upon each cycle, PCR can thereby amplify DNA fragments up to 10⁸-fold in a short period.

2) Fluorescence detection

SYBR Green fluorescence is enormously increased upon binding to double-stranded DNA. During the extension phase, more and more SYBR Green will bind to the PCR product, resulting in an increased fluorescence. Consequently, during each subsequent PCR cycle more fluorescence signal will be detected.



III. Content

Cat. No	Size	Descriptions
K-6210	96 tests	Exicycler™ 96, 8-well strip, 20 µl reaction
K-6200	96 tests	Exicycler™ 96, 8-well strip, 50 µl reaction
K-6211	96 tests	ABI 7500, 8-well strip, 20 µl reaction
K-6201	96 tests	ABI 7500, 8-well strip, 50 µl reaction
K-6212	96 tests	Opticon®, 8-well strip, 20 µl reaction
K-6202	96 tests	Opticon®, 8-well strip, 50 µl reaction
K-6213	96 tests	Exicycler™ 96, 96-well plate, 20 µl reaction
K-6203	96 tests	Exicycler™ 96, 96-well plate, 50 µl reaction
K-6214	96 tests	ABI 7500, 96-well plate, 20 µl reaction
K-6204	96 tests	ABI 7500, 96-well plate, 50 µl reaction

Cat. No	Kit Contents
K-6210	8-well strip x 12 each
K-6211	DEPC-D.W. 1.2 ml x 2 tubes
K-6212	* ROX dye (50X) 0.2 ml x 1 tube (only in K-6211)
K-6200	8-well strip x 12 each
K-6201	DEPC-D.W. 1.2 ml x 4 tubes
K-6202	* ROX dye (50X) 0.2 ml x 1 tube (only in K-6201)
K-6213	96-well plate x 1 each
K-6214	DEPC-D.W. 1.2 ml x 2 tubes
	* ROX dye (50X) 0.2 ml x 1 tube (only in K-6214)
K-6203	96-well plate x 1 each
K-6204	DEPC-D.W. 1.2 ml x 4 tubes
	* ROX dye (50X) 0.2 ml x 1 tube (only in K-6204)

* ROX dye is used for normalization of intensity by background subtraction. The use of ROX dye (50X) is recommended for Applied Biosystems 7500 Real-Time PCR System.

* The use of ROX dye is not required for Bioneer Exicycler 96 and Bio-Rad DNA engine Opticon, iCycler IQ5 real-time instruments.

IV. Storage

AccuPower GreenStar qPCR PreMix should be stored at -20°C upon received, and are stable until the expiry date stated on the label.

V. Additionally required materials & Devices

- Thermal Cycler for real-time PCR (authorized instruments)
- Target-specific primers
- Calibrated micropipette
- Sterilized micropipette tips with filters
- Optical adhesive films for real-time PCR
- High-speed Centrifuge with rotors for microtiter plates
- Vortex mixer, desktop centrifuge
- Disposable powder-free gloves

VI. General precautions

- Wear gloves during experiments to prevent contamination.
- Store positive materials, such as samples and control templates, in separated freezer from freezers for the kit.
- Add templates to the reaction mixture in clean bench or a spatially separated facility
- Place tubes / plate at room temperature at least 5 min before use.
- **Vortex and centrifuge briefly tubes before load tubes into instruments.**

AccuPower® GreenStar™ qPCR PreMix

VII. Protocol

Recommended Protocol Using *Exicycler 96* version 3.0 (Bioneer Co.), Applied Biosystems 7500 Real-time PCR System (Applied Biosystems) and DNA Engine Opticon (Formerly, MJ Research; Bio-Rad Inc.)

1. Add following PCR reagents into *AccuPower GreenStar* qPCR PreMix tube (per reaction)

	20 µl Rxn	50 µl Rxn
PCR F-Primer (10 pmole)	1-2 µl	1-2 µl
PCR R-Primer (10 pmole)	1-2 µl	1-2 µl
Template	5-10 µl	5-10 µl
DEPC-distilled water.	Adjust to 20 µl	Adjust to 50 µl

- Seal the Optical adhesive film for real-time PCR on tube or plate
- Completely mix by vigorous vortexing for resuspension of PreMix pellets.
- Centrifuge at 3,000 rpm, for 2 min
- Start Real-time PCR instrument and load it
- Program the PCR setting

Step	Condition	Cycle
Pre-Denaturation	95 °C, 1-5 min	1
Denaturation	95 °C, 5-20 sec	40-45
Annealing/Extension	55-60 °C, 40-45 sec	
Detection(Scan)		
Melting	-	1

7. After reaction is completed, perform data analysis.

*** This recommended protocol can be modified to get the optimal results, based on used real-time PCR instrument and target DNA sequences.**

VIII. Experimental Example

- Target : Envelope gene of West Nile Virus (WNV)
- Primer: Designed using Primer3 Plus & purchase from Bioneer Co. (SOUTH KOREA)
- Template: Plasmid DNA containing Envelop gene region of WNV
- Used Reagent composition (per 50 µl reaction)

WNV Forward Primer (10 pmole)	2 µl
WNV Reverse Primer (10 pmole)	2 µl
Template(10 ⁹ ~ 10 ³ copies / rxn)	5 µl
DEPC-distilled water.	Final 50 µl
- Program the PCR setting

Step	Condition	Cycle
Pre-Denaturation	95 °C, 1 min	1
Denaturation	95 °C, 5 sec	40
Annealing/Extension	55 °C, 40 sec	
Detection(Scan)		
Melting	-	1

6. Perform Data analysis




Figure 1. Data using *AccuPower GreenStar* qPCR PreMix
AccuPower GreenStar qPCR PreMix provides at least 7 orders of magnitude in dynamic range (10⁹ ~ 10³ copies/reaction).
 (A) Amplification curve, (B) Standard curve, (C) Melting curve.
 All data were obtained using *Exicycler 96* Real-time Quantitative Thermal Block (Bioneer Co.).




Figure 2. Data using various kinds of Real-time PCR Instruments

AccuPower GreenStar qPCR PreMix is applicable to most of commercially available real-time quantitative PCR instruments. West Nile Virus (WNV) primers were added into *AccuPower GreenStar* qPCR PreMix. A series of WNV positive control diluents were tested.

(A) Amplification curve, standard curve and melting curve using ABI 7500 Fast Real-time PCR machine (Applied Biosystems).

(B) Amplification curve, standard curve and melting curve using DNA Engine Opticon Real-time PCR machine (MJ Research, currently Bio-Rad Inc.).

IX. Related Products

Cat. No.	Product
K-6100 ~ K-6112	<i>AccuPower® DualStar™</i> qPCR PreMix, <i>Exicycler 96</i> , ABI 7500, Opticon 8-well strip, 96 tests /pkg
K-6103 ~ K-6114	<i>AccuPower® DualStar™</i> qPCR PreMix, <i>Exicycler 96</i> , ABI 7500 96-well plate, 96 tests /pkg
K-3032	<i>AccuPrep™</i> Genomic DNA Extraction Kit, 100 extractions
K-3033	<i>AccuPrep™</i> Viral RNA Extraction Kit, 100 extractions
A-2060	<i>Exicycler™ 96</i> Real-Time Quantitative Thermal Block