

Instruction Manual

pSV40 β Mammalian *lacZ*hls12co Expression Vector

Version 1.01
March 29, 2004

A Limited Label License covers this product (see Purchaser Notification). By use of this product, you accept the terms and conditions of the Limited Label License.

A plasmid expression vector for cloning and expression of proteins into mammalian cells with detection using the *lacZ*hls12co -Galactosidase marker gene.

Catalog no. M1018







Table of Contents

Table of Contents.....	2
Important Information.....	2
Shipping and Storage	2
Contents	2
Quality Control	2
Accessory Products	3
Methods.....	4
Using pSV40 β /acZnIs12co.....	5
Transfection.....	6
Expression and Analysis.....	7
Creating Stable Cell Lines.....	7
Appendix.....	9
Recipes	9
Map and Features of pSV40 β /acZnIs12co.....	10
Technical Service.....	18
Purchaser Notification	18
References.....	19

Important Information:

MSDS Sheets and product safety information are available by request from Marker Gene Technologies, Inc. and by accessing our web site at www.markergene.com.

Shipping and Storage

The pSV40 β /acZnIs12co vector is shipped at room temperature. Store at -20°C once resuspended. Products are guaranteed for six months from date of shipment when stored properly.

Contents

Item	Concentration
pSV40 β /acZnIs12co Vector, lyophilized in TE buffer, pH 8.0	20 μg

Quality Control

The pSV40 β /acZnIs12co vector has been qualified by restriction endonuclease digestion. pSV40 β /acZnIs12co is further qualified by transformation using an appropriate *E. coli* and mammalian cell strain in culture.



Accessory Products

Additional products that may be used with the pSV40 β -*lacZ*hls12co vector are now available from Marker Gene.

Ordering information is provided below.

Product	Unit Size	Catalog no.
pCMV β Mammalian <i>lacZ</i> Expression Vector	20 μ g	M0951
pSV40 β Mammalian <i>lacZ</i> Expression Vector	20 μ g	M0952
Expression of your recombinant fusion protein can be detected using:		
Fluorescein di- β -D-Galactopyranoside (FDG)	5 mg	M0250
Methylumbelliferyl- β -D-Galactopyranoside (MUG)	1g	M0241
Resorufin- β -D-Galactopyranoside (Res-Gal)	10 mg	M0203
Trifluoromethylumbelliferyl- β -D-Galactopyranoside	100mg	M0252
Carboxyumbelliferyl- β -D-Galactopyranoside (CUG)	5mg	M0257
FACS Fluorescent Blue <i>lacZ</i> β -Galactosidase Detection Kit	1 kit	M0255
β -Galactosidase Sample Kit	1 kit	M0276
<i>in vivo</i> <i>lacZ</i> β -Galactosidase Detection Kit	1 kit	M0259
Chemiluminescent <i>lacZ</i> β -Galactosidase Detection Kit	1 kit	M0855



Methods

Overview

Description

- This common eukaryotic expression vector, pSV40 β expresses the full-length codon-optimized β -galactosidase gene (*lacZ*nl512co) under the control of simian virus 40 (SV40) early promoter. When expressed in mammalian cells, the codon-optimized gene results in expression levels of β -galactosidase 15-fold higher than those resulting from an analogous construct containing the native *E. coli* gene sequence. Enhanced transcript stability and increased translational efficiency provide for increased β -galactosidase expression, as suggested by RNA analysis. In addition, codon-optimization results in the elimination of several cryptic splice acceptor sites that are present in the native *E. coli* gene sequence and increases the amounts of un-spliced, full-length genomic RNA when used in a lentiviral vector containing a 5' splice donor. The nl512 variant results from the addition of a twelve amino acid sequence, ProLysLysLysArgLysValGluAspProLysAsp (from the SV40 T antigen nuclear localization signal) after the methionine initiation residue.
- This vector is very useful for transfection of mammalian cells in culture and for use in other species. The β -galactosidase enzyme expression is enhanced by the SV40 late polyadenylation signal.
- pSV40 β /*lacZ*nl512co expression vector also contains the β -lactamase gene, which acts as a selection marker (100 μ g/mL ampicillin resistance) in *E. coli* host.
- The β -galactosidase gene can be excised using the 5' *Xho*I and 3' *Not*I sites to allow the insertion of other genes to be expressed under the same regulatory elements in mammalian cells.
- For a map of pSV40 β , see page 10.

The pSV40 β /*lacZ*nl512co Vector System

The pSV40 β /*lacZ*nl512co vector is a cloning vector that takes advantage of the site-specific recombination properties of bacteriophage lambda (Landy, 1989) to provide a rapid and highly efficient way to move your gene of interest into multiple vector systems. To express your gene of interest simply:

1. Clone your gene of interest into an entry vector to create an entry clone.
2. Perform an LR recombination reaction between the entry clone and a destination vector (e.g. pSV40 β) to generate the expression vector.
3. Transfect your expression clone into the cell line of choice for stable expression of your gene of interest.



Using pSV40 β /acZnIs12co

The pSV40 β /acZnIs12co vector is supplied as a supercoiled plasmid. Although Marker Gene has previously recommended using a linearized destination vector for more efficient recombination, further testing has found that linearization of this vector is NOT required to obtain optimal results for a downstream application.

Propagating pSV40 β :

If you wish to propagate and maintain pSV40 β /acZnIs12co, we recommend using JM109 *E. coli* Competent Cells.

Entry Clone:

To recombine your gene into pSV40 β /acZnIs12co, you should have an entry clone containing your gene of interest.

Points to consider before recombining:

Your insert should contain a Kozak consensus sequence with an ATG initiation codon for proper initiation of translation (Kozak, 1987; Kozak, 1991; Kozak, 1990). Other sequences are also possible, but the G or A at position -3 and the G at position +4 are the most critical for functional expression. If you wish to include the V5 epitope and 6xHis tag, your gene in the entry clone should not contain a stop codon. The gene should also be designed to be in frame with the C-terminal epitope tag after recombination. If you do NOT wish to include the V5 epitope and 6xHis tag, please be sure that your gene contains a stop codon in the entry clone.



Transfection

Introduction: This section provides general information for transfecting your expression clone into the mammalian cell line of choice. We recommend that you include a positive control vector and a mock transfection (negative control) in your experiments to evaluate your results.

Plasmid Preparation:

Once you have generated your expression clone, you must isolate plasmid DNA for transfection. Plasmid DNA must be very clean and free from phenol and sodium chloride. Contaminants will kill the cells, and salt will interfere with lipid complexing, decreasing transfection efficiency. We recommend isolating the plasmid using a mini/midi prep kit such as the one available from Qiagen (Plasmid Midi Kit, Cat# 12143).

Methods of Transfection:

For established cell lines (e.g. HeLa), consult original references or the supplier of your cell line for the optimal method of transfection. We recommend that you follow exactly the protocol for your particular cell line. Please pay particular attention to cell medium requirements, confluency and when to pass the cells, and at what dilution to split the cells. Further information is provided in *Current Protocols in Molecular Biology* (Ausubel *et al.*, 1994). Methods for transfection include calcium phosphate (Chen and Okayama, 1987; Wigler *et al.*, 1977), lipid-mediated (Felgner *et al.*, 1989; Felgner and Ringold, 1989) and electroporation (Chu *et al.*, 1987; Shigekawa and Dower, 1988). If you wish to use a lipid-based reagent for transfection, we recommend using Lipofectamine Reagent available from GIBCO-BRL. For more information contact our Technical Assistance Staff (www.markergene.com or techservice@markergene.com).

Positive Control:

We recommend the use of a positive control vector for mammalian cell transfection and expression which may be used to optimize recombinant protein expression levels in your particular cell line. A vector that allows expression of a C-terminally tagged β -galactosidase fusion protein that may be detected by Western blot or functional assay provides the easiest way to measure protein expression levels. Consult our technical assistance for more information about C-terminal fusion protein expression systems.

To propagate and maintain the plasmid:

1. Resuspend the vector in 20 μ l sterile water to prepare a 1 μ g/ μ L stock solution and store at -20°C. Use the stock solution to transform a *recA*, *endA* *E. coli* strain like TOP10, JM109, or equivalent.
2. Select transformants on LB agar plates containing 50-100 μ g/ml ampicillin.
3. Prepare a glycerol stock of a transformant containing plasmid for long-term storage.



Expression and Analysis

Expression of your gene of interest from the expression clone can be performed in transiently transfected cells or stable cell lines.

Preparation of Cell Lysates:

To lyse cells:

1. Wash cell monolayers ($\sim 5 \times 10^5$ to 1×10^6 cells) once with phosphate-buffered saline (PBS, available from Gibco, Catalog no. 10010-023 or see page 9 for a recipe).
2. Scrape cells into 1 ml PBS and pellet the cells at $1500 \times g$ for 5 minutes.
3. Resuspend in 50 μ l Cell Lysis Buffer (Marker Gene Product # M0626-003 or see page 9 for a recipe). Other cell lysis buffers are also suitable. Vortex mix.
4. Incubate cell suspension at 37°C for 10 minutes to lyse the cells. Note: You may prefer to lyse the cells at room temperature or on ice if degradation of your protein is a potential problem.
5. Centrifuge the cell lysate at $10,000 \times g$ for 10 minutes at 4°C to pellet nuclei and transfer the supernatant to a fresh tube. Assay the lysate for protein concentration. Do not use protein assays utilizing Coomassie Blue or other dyes. NP-40 may interfere with the binding of the dye with the protein.
6. Add SDS-PAGE sample buffer (see page 12 for a recipe) to a final concentration of 1X and boil the sample for 5 minutes.
7. Load 20 μ g of lysate onto an SDS-PAGE gel and electrophorese. Use the appropriate percentage of acrylamide to resolve your fusion protein.

Detecting Recombinant Fusion Proteins:

To detect expression of your recombinant fusion protein by Western blot analysis, you may use the Anti-V5 antibodies or the Anti-His(C-term) antibodies available from Invitrogen or Amersham Biosciences or an antibody to your protein of interest.

Assay for β -galactosidase:

If you use a positive control vector, you may assay for β -galactosidase expression by Western blot analysis or activity assay using cell-free lysates (Miller, 1972). X-Gal staining, or fluorescence detection are common methods of analysis. Marker Gene offers a FACS Fluorescent Blue *lacZ* β -Galactosidase Detection Kit (Product M0255), and the *in vivo lacZ* β -Galactosidase Detection Kit (Product M0259) for fast and easy detection of β -galactosidase expression.

Purification of Recombinant Fusion Proteins:

The presence of the C-terminal polyhistidine (6xHis) tag in your recombinant fusion protein allows for purification using a metal-chelating resin (available from Life Technologies). Note: Other purification methods are suitable.

Creating Stable Cell Lines:

The neomycin resistance gene can be cloned into the pSV40 β *lacZ*hls12co vector to allow selection of stable cell lines. If you wish to create stable cell lines, transfect your construct into the mammalian cell line of choice and select for foci using Geneticin®. General guidelines are provided below.

To obtain stable transfectants, we recommend that you linearize your pSV40 β *lacZ*hls12co construct before transfection. While linearizing the vector may not improve the efficiency of transfection, it increases the chances that the vector does not integrate in a way that disrupts elements necessary for expression in mammalian cells. To linearize your construct, cut at a unique site that is neither located within a critical element nor within your gene of interest.



Geneticin® (G418) blocks protein synthesis in mammalian cells by interfering with ribosomal function. It is an aminoglycoside, similar in structure to neomycin, gentamycin, and kanamycin. Expression in mammalian cells of the bacterial aminoglycoside phosphotransferase gene (APH), derived from *Tn5*, results in detoxification of Geneticin® (Southern and Berg, 1982).

Geneticin® Selection Guidelines:

Geneticin® is available from GIBCO (Catalog no. 11811-023). Use as follows:

1. Prepare Geneticin® in a buffered solution (e.g. 100 mM HEPES, pH 7.3).
2. Use 100 to 1000 µg/ml of Geneticin® in complete medium.
3. Calculate concentration based on the amount of active drug.
4. Test varying concentrations of Geneticin® on your cell line to determine the concentration that kills your cells (kill curve). Cells differ in their susceptibility to Geneticin®. Cells will divide once or twice in the presence of lethal doses of Geneticin® Selective Antibiotic, so the effects of the drug takes several days to become apparent. Complete selection for positive clones of cells can take up to 2 to 3 weeks of growth in selection medium.



Appendix

Recipes

LB (Luria-Bertani) Medium and Plates Composition:

1.0% Tryptone

0.5% Yeast Extract

1.0% NaCl

pH 7.0

1. For 1 liter, dissolve 10 g tryptone, 5 g yeast extract, and 10 g NaCl in 950 ml deionized water.

2. Adjust the pH of the solution to 7.0 with NaOH and bring the volume up to 1 liter.

3. Autoclave on liquid cycle for 20 minutes at 15 psi. Allow solution to cool to 55°C and add antibiotic if needed. Store at room temperature or at +4°C.

LB agar plates:

1. Prepare LB medium as above, but add 15 g/L agar before autoclaving.

2. Autoclave on liquid cycle for 20 minutes at 15 psi.

3. After autoclaving, cool to ~55°C, add antibiotic if needed, and pour into 10 cm plates.

4. Let harden, then invert and store at +4°C.

Cell Lysis Buffer:

25 mM Tris-phosphate (pH 7.8) containing 10% glycerol, 1% Triton X-100, 1 mg/ml BSA, 2 mM EGTA and 2 mM DTT 50 mM Tris, pH 7.8.

4X SDS-PAGE Sample Buffer:

1. Combine the following reagents:

0.5 M Tris-HCl, pH 6.8, 5 ml

Glycerol (100%), 4 ml

β -mercaptoethanol, 0.8 ml

Bromophenol Blue, 0.04 g

SDS, 0.8 g


2. Bring the volume to 10 ml with sterile water.

3. Aliquot and freeze at -20°C until needed.



Map and Features of pSV40 β /lacZnls12co:

The map below shows the elements of pSV40 β /lacZnls12co. The β -galactosidase gene can be excised using the 5' *Xho*I and 3' *Not*I sites to allow the insertion of other genes to be expressed under the same regulatory elements in mammalian cells.



Full Length Sequence Data of pSV40 β /lacZnls12co:

```

1  GAATTCTAGT  TGTGGTTTGT  CCAAATCAT  CAATGTATCT  TATCATGTCT
51  GGATCCGCTG  TGGAATGTGT  GTCAGTTAGG  GTGTGGAAAG  TCCCCAGGCT
101  CCCCAGCAGG  CAGAAGTATG  CAAAGCATGC  ATCTCAATTA  GTCAGCAACC
151  AGGTGTGGAA  AGTCCCCAGG  CTCGCCAGCA  GGCAGAAGTA  TGCAAAGCAT
201  GCATCTCAAT  TAGTCAGCAA  CCATAGTCCC  GCCCCTAACT  CCGCCCATCC
251  CGCCCCTAAC  TCCGCCAGT  TCCGCCATT  CTCCGCCCA  TGGCTGACTA
301  ATTTTTTTTT  TTTATGCAGA  GGCCGAGGCC  GCCTCGGCCT  CTGAGCTATT
351  CCAGAAGTAG  TGAGGAGGCT  TTTTGGAGG  CCTAGGCTTT  TGCAAAAAGC
401  TTGGACACAA  GACAGGCTTG  CGAGATATGT  TTGAGAATAC  CACTTTATCC
451  CGCGTCAGGG  AGAGGCAGTG  CGTAAAAAGA  CGCGGACTCA  TGTGAAATAC
501  TGGTTTTTAG  TGCGCCAGAT  CTCGAGGTCG  ACGGTATCGA  TAAGCTTAAC

```



551 CATCAGCAAG CAGGTCATTG TGCCACCACC ATGCCCAAGA AGAAGAGGAA
601 GGTGGAGGAC CCCAAGGACA TCACCGACTC CCTGGCCGTG GTGCTGCAGC
651 GCCGCGACTG GGAGAACCCC GCGGTGACCC AGCTGAACCG CCTGGCCGCC
701 CACCCCCCCT TCGCCTCCTG GCGCAACTCC GAGGAGGCCG GCACCGACCG
751 CCCCTCCAG CAGCTGCGCT CCCTGAACGG CGAGTGGCGC TTCGCCTGGT
801 TCCCCGCCCC CGAGGCCGTG CCCGAGTCCT GGCTGGAGTG CGACCTGCCC
851 GAGGCCGACA CCGTGGTGGT GCCCTCCAAC TGGCAGATGC ACGGCTACGA
901 CGCCCCCATC TACACCAACG TGACCTACCC CATCACCGTG AACCCCCCTT
951 TCGTGCCAC CGAGAACCCC ACCGGCTGCT ACTCCCTGAC CTTCAACGTG
1001 GACGAGTCCT GGCTGCAGGA GGGCCAGACC CGCATCATCT TCGACGGCGT
1051 GAACTCCGCC TTCCACCTGT GGTGCAACGG CCGCTGGGTG GGCTACGGCC
1101 AGGACTCCCG CCTGCCCTCC GAGTTCGACC TGTCCGCCTT CCTGCGCGCC
1151 GGCGAGAACC GCCTGGCCGT GATGGTGCTG CGCTGGTCCG ACGGCTCCTA
1201 CCTGGAGGAC CAGGACATGT GGCGCATGTC CGGCATCTTC CGCGACGTGT
1251 CCCTGCTGCA CAAGCCCACC ACCCAGATCT CCGACTTCCA CGTGGCCACC
1301 CGCTTCAACG ACGACTTCTC CCGCGCCGTG CTGGAGGCCG AGGTGCAGAT
1351 GTGCGGCGAG CTGCGCGACT ACCTGCGCGT GACCGTGTCC CTGTGGCAGG
1401 GCGAGACCCA GGTGGCCTCC GGCACCGCCC CCTTCGGCGG CGAGATCATC
1451 GACGAGCGCG GCGGCTACGC CGACCGCGTG ACCCTGCGCC TGAACGTGGA
1501 GAACCCCAAG CTGTGGTCCG CCGAGATCCC CAACCTGTAC CGCGCCGTGG
1551 TGGAGCTGCA CACCGCCGAC GGCACCCTGA TCGAGGCCGA GGCTGCGAC
1601 GTGGGCTTCC GCGAGGTGCG CATCGAGAAC GGCTGTCTGC TGCTGAACGG
1651 CAAGCCCCTG CTGATCCGCG CGGTGAACCG CCACGAGCAC CACCCCCTGC
1701 ACGGCCAGGT GATGGACGAG CAGACCATGG TGCAGGACAT CCTGCTGATG
1751 AAGCAGAACA ACTTCAACGC CGTGCCTGTC TCCCACTACC CCAACCACCC
1801 CCTGTGGTAC ACCCTGTGCG ACCGCTACGG CCTGTACGTG GTGGACGAGG
1851 CCAACATCGA GACCCACGGC ATGGTGCCCA TGAACGCCT GACCGACGAC
1901 CCCCCTGGC TGCCCGCCAT GTCCGAGCGC GTGACCCGCA TGGTGCAGCG
1951 CGACCGCAAC CACCCCTCCG TGATCATCTG GTCCCTGGGC AACGAGTCCG
2001 GCCACGGCGC CAACCACGAC GCCCTGTACC GCTGGATCAA GTCCGTGGAC
2051 CCCTCCCGCC CCGTGCAGTA CGAGGGCGGC GGCGCCGACA CCACCGCCAC
2101 CGACATCATC TGCCCCATGT ACGCCCGCGT GGACGAGGAC CAGCCCTTCC
2151 CCGCCGTGCC CAAGTGGTCC ATCAAGAAGT GGCTGTCCCT GCCC GGCGAG
2201 ACCCGCCCCC TGATCCTGTG CGAGTACGCC CACGCCATGG GCAACTCCCT
2251 GGGCGGCTTC GCCAAGTACT GGCAGGCCTT CCGCCAGTAC CCCCCTGC
2301 AGGGCGGCTT CGTGTGGGAC TGGGTGGACC AGTCCCTGAT CAAGTACGAC
2351 GAGAACGGCA ACCCCTGGTC CGCCTACGGC GGCAGCTTCG GCGACACCCC
2401 CAACGACCGC CAGTTCTGCA TGAACGGCCT GGTGTTCCGCC GACCGCACCC
2451 CCCACCCCGC CCTGACCGAG GCCAAGCACC AGCAGCAGTT CTTCCAGTTC
2501 CGCCTGTCCG GCCAGACCAT CGAGGTGACC TCCGAGTACC TGTTCCGCCA
2551 CTCCGACAAC GAGCTGCTGC ACTGGATGGT GGCCCTGGAC GGCAAGCCCC
2601 TGGCCTCCGG CGAGGTGCCC CTGGACGTGG CCCCCAGGG CAAGCAGCTG
2651 ATCGAGCTGC CCGAGCTGCC CCAGCCCGAG TCCGCCGGCC AGCTGTGGCT
2701 GACCGTGC GC GTGGTGCAGC CCAACGCCAC CGCCTGGTCC GAGGCCGGCC
2751 ACATCTCCGC CTGGCAGCAG TGGCGCCTGG CCGAGAACCT GTCCGTGACC
2801 CTGCCC GCCG CCTCCCACGC CATCCCCAC CTGACCACCT CCGAGATGGA
2851 CTTCTGCATC GAGCTGGGCA ACAAGCGCTG GCAGTTCAAC CGCCAGTCCG
2901 GCTTCTGTG CCAGATGTGG ATCGGCGACA AGAAGCAGCT GCTGACCCCC
2951 CTGCGCGACC AGTTCACCCG CGCCCCCTG GACAACGACA TCGGCGTGTG
3001 CGAGGCCACC CGCATCGACC CCAACGCCTG GGTGGAGCGC TGGAAGGCCG
3051 CCGGCCACTA CCAGGCCGAG GCCGCCCTGC TGCAGTGCAC CGCCGACACC
3101 CTGGCCGACG CCGTGTCTGAT CACCACCGCC CACGCCTGGC AGCACCAGGG
3151 CAAGACCCTG TTCATCTCCC GCAAGACCTA CCGCATCGAC GGCTCCGGCC
3201 AGATGGCCAT CACCGTGGAC GTGGAGGTGG CCTCCGACAC CCCCCACCCC
3251 GCCCGCATCG GCCTGAACTG CCAGCTGGCC CAGGTGGCCG AGCGCGTGAA



3301 CTGGCTGGGC CTGGGCCCCC AGGAGAACTA CCCCAGCCGC CTGACCGCCG
3351 CCTGCTTCGA CCGCTGGGAC CTGCCCCTGT CCGACATGTA CACCCCCTAC
3401 GTGTTCCCTT CCGAGAACGG CCTGCGCTGC GGCACCCGCG AGCTGAACTA
3451 CGGCCCCAC CAGTGGCGCG GCGACTTCCA GTTCAACATC TCCCGCTACT
3501 CCCAGCAGCA GCTGATGGAG ACCTCCCACC GCCACCTGCT GCACGCCGAG
3551 GAGGGCACCT GGCTGAACAT CGACGGCTTC CACATGGGCA TCGGCGGCGA
3601 CGACTCCTGG TCCCCCTCCG TGTCCGCCGA GTTCCAGCTG TCCGCCGGCC
3651 GCTACCACTA CCAGCTGGTG TGGTGCCAGA AGTAGTCTAG AGCGGCCGCG
3701 GGGATCCAGA CATGATAAGA TACATTGATG AGTTTGGACA AACCACAACCT
3751 AGAATGCAGT GAAAAAATG CTTTATTTGT GAAATTTGTG ATGCTATTGC
3801 TTTATTTGTA ACCATTATAA GCTGCAATAA ACAAGTTAAC AACAACAATT
3851 GCATTCATTT TATGTTTCAG GTTCAGGGGG AGGTGTGGGA GGTTTTTTTCG
3901 GATCCTCTAG AGTCGACCTG CAGGCATGCA AGCTTGGCGT AATCATGGTC
3951 ATAGCTGTTT CCTGTGTGAA ATTGTTATCC GCTCACAATT CCACACAACA
4001 TACGAGCCGG AAGCATAAAG TGTAAGCCTT GGGGTGCCTA ATGAGTGAGC
4051 TAACTCACAT TAATTGCGTT GCGCTCACTG CCCGCTTTCC AGTCGGGAAA
4101 CCTGTCTGTC CAGCTGCATT AATGAATCGG CCAACGCGCG GGGAGAGGCG
4151 GTTTGCGTAT TGGGCGCTCT TCCGCTTCCT CGCTCACTGA CTCGCTGCGC
4201 TCGGTCGTTT GGCTGCGGCG AGCGGTATCA GCTCACTCAA AGGCGGTAAT
4251 ACGGTTATCC ACAGAATCAG GGGATAACGC AGGAAAGAAC ATGTGAGCAA
4301 AAGGCCAGCA AAAGGCCAGG AACCCTAAAA AGGCCGCGTT GCTGGCGTTT
4351 TTCCATAGGC TCCGCCCCC TGACGAGCAT CACAAAAATC GACGCTCAAG
4401 TCAGAGTTGG CGAAACCCGA CAGGACTATA AAGATACCAG GCGTTTCCCC
4451 CTGGAAGCTC CCTCGTGCGC TCTCCTGTTT CGACCCTGCC GCTTACCGGA
4501 TACCTGTCCG CCTTTCTCCC TTCGGGAAGC GTGGCGCTTT CTCATAGCTC
4551 ACGCTGTAGG TATCTCAGTT CGGTGTAGGT CGTTCGCTCC AAGCTGGGCT
4601 GTGTGCACGA ACCCCCCGTT CAGCCCGACC GCTGCGCCTT ATCCGGTAAC
4651 TATCGTCTTG AGTCCAACCC GGTAAGACAC GACTTATCGC CACTGGCAGC
4701 AGCCACTGGT AACAGGATTA GCAGAGCGAG GTATGTAGGC GGTGCTACAG
4751 AGTTCTTGAA GTGGTGGCCT AACTACGGCT AACTAGAAG GACAGTATTT
4801 GGTATCTGCG CTCTGCTGAA GCCAGTTACC TTCGGAAGAA GAGTTGGTAG
4851 CTCTTGATCC GGCAACAAA CCACCGCTGG TAGCGGTGGT TTTTTTGTTC
4901 GCAAGCAGCA GATTACGCGC AGAAAAAAG GATCTCAAGA AGATCCTTTG
4951 ATCTTTTCTA CGGGGTCTGA CGCTCAGTGG AACGAAAAC CACGTTAAGG
5001 GATTTTGGTC ATGAGATTAT CAAAAAGGAT CTTACCTAG ATCCTTTTAA
5051 ATTAAAAATG AAGTTTTTAA TCAATCTAAA GTATATATGA GTAAACTTGG
5101 TCTGACAGTT ACCAATGCTT AATCAGTGAG GCACCTATCT CAGCGATCTG
5151 TCTATTTCTG TCATCCATAG TTGCCTGACT CCCCCTCGTG TAGATAACTA
5201 CGATACGGGA GGGCTTACCA TCTGGCCCCA GTGCTGCAAT GATACCGCGA
5251 GACCCACGCT CACCGGCTCC AGATTTATCA GCAATAAAC AGCCAGCCGG
5301 AAGGCCCGAG CGCAGAAGTG GTCTGCAAC TTTATCCGCC TCCATCCAGT
5351 CTATTAATTG TTGCCGGGAA GCTAGAGTAA GTAGTTCGCC AGTTAATAGT
5401 TTGCGCAACG TTGTTGCCAT TGCTACAGGC ATCGTGGTGT CACGCTCGTC
5451 GTTTGGTATG GCTTCATTCA GCTCCGGTTC CCAACGATCA AGGCGAGTTA
5501 CATGATCCCC CATGTTGTGC AAAAAAGCGG TTAGCTCCTT CGGTCTCCG
5551 ATCGTTGTCA GAAGTAAGTT GGCCGAGTG TTATCACTCA TGGTTATGGC
5601 AGCACTGCAT AATTCTCTTA CTGTCATGCC ATCCGTAAGA TGCTTTTCTG
5651 TGACTGGTGA GACTCAACC AAGTCATTCT GAGAATAGTG TATGCGGCGA
5701 CCGAGTTGCT CTTGCCCGGC GTCAATACGG GATAATACCG CGCCACATAG
5751 CAGAACTTTA AAAGTGCTCA TCATTGAAA ACGTTCTTCG GGGCGAAAAC
5801 TCTCAAGGAT CTTACCGCTG TTGAGATCCA GTTCGATGTA ACCCACTCGT
5851 GCACCCAAC GATCTTCAGC ATCTTTTACT TTCACCAGCG TTTCTGGGTG
5901 AGCAAAAACA GGAAGGCAA ATGCCGCAA AAAGGGAATA AGGGCGACAC
5951 GGAAATGTTG AATACTCATA CTCTTCCTTT TTCAATATTA TTGAAGCATT
6001 TATCAGGGTT ATTGTCTCAT GAGCGGATAC ATATTTGAAT GTATTTAGAA



```
6051 AAATAAACAA ATAGGGGTTC CGCGCACATT TCCCCGAAAA GTGCCACCTG
6101 ACGTCTAAGA AACCATATT ATCATGACAT TAACCTATAA AAATAGGCGT
6151 ATCACGAGGC CCTTTCGTCT CGCGCGTTTC GGTGATGACG GTGAAAACCT
6201 CTGACACATG CAGCTCCCGG AGACGGTCAC AGCTTGTCTG TAAGCGGATG
6251 CCGGGAGCAG ACAAGCCCGT CAGGGCGCGT CAGCGGGTGT TGGCGGGTGT
6301 CGGGGCTGGC TTAACATATG GGCATCAGAG CAGATTGTAC TGAGAGTGCA
6351 CCATATGCGG TGTGAAATAC CGCACAGATG CGTAAGGAGA AAATACCGCA
6401 TCAGGCGCCA TTCGCCATTC AGGCTGCGCA ACTGTTGGGA AGGGCGATCG
6451 GTGCGGGCCT CTTCGCTATT ACGCCAGCTG GCGAAAGGGG GATGTGCTGC
6501 AAGGCGATTA AGTTGGGTAA CGCCAGGGTT TTCCAGTCA CGACGTTGTA
6551 AAACGACGGC CAGT
```

Restriction Sites:

AarI (CACCTGCnnnn'nnnn_)

Cuts 1 time.

Cuts at position 3543.

AatII (G_ACGT'C) [ZraI]

Cuts 1 time.

Cuts at position 6104.

AccI (GT'mk_AC) [FblI,XmiI]

Cuts 2 times.

AcII (AA'CG_TT) [Psp1406I]

Cuts 2 times.

AcuI (CTGAAGnnnnnnnnnnnnnn'_nn') [Eco57I]

Cuts 2 times.

AfeI (AGC'GCT) [Eco47III,Aor51HI,FunI]

Cuts 2 times.

AhdI (GACnn_n'nnGTC) [Eam1105I,AspEI,DriI,EclHKI]

Cuts 2 times.

AleI (CACnn'nnGTG) [OliI]

Cuts 2 times.

ApaI (G_GGCC'C) [Bsp120I,PspOMI]

[dcm methylated]

Cuts 1 time.

Cuts at position [3317].

ApoI (r'AATT_y) [AcsI,XapI]

Cuts 2 times.



AvrII (C'CTAG_G) [AspA2I,BlnI,XmaJI]

Cuts 1 time.

Cuts at position 381.

BanII (G_rGCy'C) [Eco24I,EcoT38I,FriOI]

Cuts 1 time.

Cuts at position 3317.

BciVI (GTATCCnnnnn_n') [BfuI]

Cuts 2 times.

BfrBI (ATG'CAT) [EcoT22I,Mph1103I,NsiI,Zsp2I]

Cuts 2 times.

BglII (A'GATC_T)

Cuts 2 times.

BsaBI (GATnn'nnATC) [Bse8I,BseJI,MamI]

[dam methylated]

Cuts 2 times.

BsmI (GAATG_Cn') [BsaMI,Mva1269I,PctI]

Cuts 2 times.

BsmBI (CGTCTCn'nnnn_) [Esp3I]

Cuts 2 times.

BsrDI (GCAATG_nn') [Bse3DI,BseMI]

Cuts 2 times.

BsrGI (T'GTAC_A) [Bsp1407I,BstAUI,SspBI]

Cuts 1 time.

Cuts at position 3387.

BssHII (G'CGCG_C) [BsePI,PauI]

Cuts 1 time.

Cuts at position 1144.

BstEII (G'GTnAC_C) [BstPI,Eco91I,EcoO65I,PspEI]

Cuts 1 time.

Cuts at position 2524.

BstXI (CCAn_nnnn'nTGG)

Cuts 1 time.

Cuts at position 3214.



ClaI (AT'CG_AT)
[BanIII, Bsa29I, BseCI, Bsp106I, BspDI, BspXI, Bsu15I, BsuTUI, ZhoI]
[dam methylated]
Cuts 1 time.
Cuts at position 537.

DraIII (CAC_nnn'GTG) [AdeI]
Cuts 1 time.
Cuts at position 1815.

DrdI (GACnn_nn'nnGTC) [AasI, DseDI]
Cuts 2 times.

EcoRI (G'AATT_C) [FunII]
Cuts 1 time.
Cuts at position 1.

FalI (AAGnnnnnnCTTnnnnnnnnn_nnnnn')
Cuts 1 time.
Cuts at position 432.

FseI (GG_CCGG'CC)
Cuts 1 time.
Cuts at position 2748.

FspAI (rTGC'GCAy)
Cuts 1 time.
Cuts at position 1619.

HpaI (GTT'AAC) [KspAI]
Cuts 1 time.
Cuts at position 3837.

MfeI (C'AATT_G) [MunI]
Cuts 1 time.
Cuts at position 3846.

MscI (TGG'CCA) [BalI, MlsI, MluNI, Msp20I]
[dcm methylated]
Cuts 2 times.

NdeI (CA'TA_TG) [FauNDI]
Cuts 1 time.
Cuts at position 6353.

NotI (GC'GGCC_GC) [CciNI]
Cuts 1 time.
Cuts at position 3693.



NsiI (A_TGCA'T) [BfrBI,EcoT22I,Mph1103I,Zsp2I]
Cuts 2 times.

PfoI (T'CCnGG_A)
Cuts 1 time.
Cuts at position 6215.

PmlI (CAC'GTG) [PmaCI,AcvI,BbrPI,Eco72I,PspCI]
Cuts 1 time.
Cuts at position 1291.

PpuMI (rG'GwC_Cy) [PpuXI,Psp5II,PspPPI]
[dcm methylated]
Cuts 2 times.

PsiI (TTA'TAA)
Cuts 1 time.
Cuts at position 3817.

PspOMI (G'GGCC_C) [ApaI,Bsp120I]
Cuts 1 time.
Cuts at position 3313.

PvuI (CG_AT'CG) [BspCI,Ple19I]
Cuts 2 times.

SacII (CC_GC'GG) [Cfr42I,KspI,Sfr303I,SgrBI]
Cuts 2 times.

SalI (G'TCGA_C)
Cuts 2 times.

SapI (GCTCTTCn'nnn_)
Cuts 1 time.
Cuts at position 4173.

SbfI (CC_TGCA'GG) [Sse8387I,SdaI]
Cuts 2 times.

ScaI (AGT'ACT) [ZrmI]
Cuts 2 times.

SexAI (A'CCwGG_T) [MabI]
[dcm methylated]
Cuts 1 time.
Cuts at position [148].



SfiI (GGCCn_nnn'nGGCC)
[dcm methylated]
Cuts 2 times.

SgrAI (Cr'CCGG_yG)
Cuts 1 time.
Cuts at position 1148.

SspI (AAT'ATT)
Cuts 1 time.
Cuts at position 5986.

TfiI (G'AwT_C) [PfeI]
Cuts 2 times.

Tth111I (GACn'n_nGTC) [AspI,PflFI,PsyI,TelI]
Cuts 1 time.
Cuts at position 1384.

XbaI (T'CTAG_A)
[dam methylated]
Cuts 2 times.

XcmI (CCAnnnn_n'nnnnTGG)
Cuts 2 times.

XhoI (C'TCGA_G) [BssHI,Paer7I,Sfr274I,SlaI,TliI]
Cuts 1 time.
Cuts at position 521.

XmnI (GAAnn'nnTTC) [Asp700I,MroXI,PdmI]
Cuts 2 times.

ZraI (GAC'GTC) [AatII]
Cuts 1 time.
Cuts at position 6102.



Technical Service

Please visit our website at www.markergene.com
Or write to techservice@markergene.com

Marker Gene is committed to providing our customers with premium-quality goods and services. We want to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about our products or service, please contact us at:

Marker Gene Technologies, Inc.
1850 Millrace Dr.
Eugene, OR 97403-1992
Toll Free: 1-888-218-4062
Local: 1-541-342-3760
Fax: 1-541-342-1960

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). No license is conveyed under the foregoing patents to use this product with any recombination sites other than those purchased from Marker Gene Technologies, Inc. or its authorized distributor. The buyer cannot modify the recombination sequence(s) contained in this product for any purpose. The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for commercial purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for the commercial purposes of the buyer, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for commercial purposes. Commercial purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. The user warrants to Marker Gene Technologies, Inc. that it will not use the purchase products or any part thereof for any purpose that is in contravention of any applicable law or regulation or in a manner that will infringe the copyright, trademarks, service marks or other intellectual property rights of third parties or violate the privacy, publicity or other personal rights of others or in any defamatory, obscene, threatening, abusive or hateful manner. The user shall indemnify Marker Gene Technologies, Inc., its employees and agents, harmless against all claims, liability, losses, damages and expenses including, without limitation, legal fees and costs arising out of or incurred as a result of any claims made, or litigation brought, against Marker Gene Technologies, Inc., its employees and agents, as a result of the use by the user of the purchased products or any part thereof for whatever purpose, to the fullest extent permitted by law.



References

- Anson DS, Limberis M. (2004) "An improved β -galactosidase reporter gene." J of Biotech 108: 17-30.
- Ausubel FM, Brent R, Kingston RE, Moore DD, Seidman JG, Smith JA, Struhl K. (1994). Current Protocols in Molecular Biology (New York: Greene Publishing Associates and Wiley-Interscience).
- Chen C, Okayama H. (1987). "High-efficiency transformation of mammalian cells by plasmid DNA." Mol Cell Biol **7**: 2745-2752.
- Chu G, Hayakawa H, Berg P. (1987). "Electroporation for the efficient transfection of mammalian cells with DNA." Nucleic Acids Res **15**: 1311-1326.
- Felgner PL, Ringold GM. (1989). "Cationic liposome-mediated transfection." Nature **337**: 387-388.
- Landy A. (1989). "Dynamic, structural, and regulatory aspects of lambda site-specific recombination." Annu Rev Biochem **58**: 913-949.
- MacGregor GR, Caskey CT. (1989). "Construction of plasmids that express *E.coli* beta-galactosidase in mammalian cells." Nucleic Acids Res **17**: 2365.
- Miller JH. (1972). Experiments in Molecular Genetics (Cold Spring Harbor, New York: Cold Spring Harbor Laboratory).
- Nolan GP, Fiering S, Nicolas JF, Herzenberg LA. (1988). "Fluorescence-activated cell analysis and sorting of viable mammalian cells based on β -D-galactosidase activity after transduction of *Escherichia coli* lacZ." Proc. Natl. Acad. Sci (USA). **85**(8):2603-2607.
- Norton PA, Coffin JM. (1984). "Bacterial β -galactosidase as a marker of Rous sarcoma virus gene expression and replication." Mol Cell Biol **5**: 281.
- Southern PJ, Berg P. (1982). "Transformation of mammalian cells to antibiotic resistance with a bacterial gene under control of the SV40 Early Region Promoter." J Molec Appl Gen **1**: 327-339.

Marker Gene Technologies, Inc.
University of Oregon Riverfront Research Park
1850 Millrace Drive
Eugene, Oregon 97403
1-888-218-4062
Email: techservice@markergene.com