



DATA SHEET

High Sensitive MiRNA Northern Blot Assay Kit

Catalog Number NB-1001

(For Research Use Only)

Introduction

Abundance of mammalian miRNAs varies dramatically in tissues. For example, expression of miR-133 is high in skeletal muscle, moderate in heart, and no detection in other tissues. To detect low expression of miRNAs requires high sensitive methods. Signosis' high sensitive miRNA northern blot assay is 10 to 100 times more sensitive than conventional biotin-based chemiluminescent detection.

Principle

RNA samples are separated through gel electrophoresis and transferred onto a membrane. Expression of a specific miRNA is detected with a biotin-labeled probe, containing two moieties – complementary sequence of the miRNA and a tag sequence. The tag sequence is then detected by an amplifier enriched with biotin molecules.

Materials provided with the kit

1. 30 μ l Ready-to-use small sizes of molecular standards (20 nt & 60 nt) (-20°C)
2. 90 μ l Gel loading buffer (RT)
3. Two 15% TBE Urea-gel (4°C)
4. Two membranes (RT)
5. 35ml 1x NB Hybridization buffer (2, RT)
6. 40ml 5x NB Hybridization wash buffer (2, RT)
7. 60ml Blocking buffer (RT)
8. 60 μ l Streptavidin-HRP conjugate (4°C)
9. 50ml 5x Detection wash buffer (RT)
10. 3.6 ml Substrate A (4°C)
11. 3.6 ml Substrate B (4°C)

Materials and equipment are needed

1. Biotin labeled miRNA probe and amplifier (HP-0XXX)
2. TBE
3. Bio-Rad gel apparatus
4. Power supplies
5. Stratagene UV cross-linker
6. Hybridization oven
7. Hybridization tubes
8. Shaker
9. Imaging system or X-ray film

Reagent preparation before experiment

- Dilute the 5x Hybridization wash buffer and 5x Detection washing buffer to 1x buffer
- 1x Hybridization wash buffer:
30ml 5x Hybridization wash buffer
120ml ddH₂O
- 1x Detection wash buffer:
40ml 5x Detection wash buffer
160 ml ddH₂O
- Prewarm 1x NB Hybridization buffer, 1x Hybridization wash buffer at 42°C for 2 hour or until the buffers are clear without visible precipitation before using.

Assay Procedure

1. Gel electrophoresis

- (1) Remove the comb and the bottom plastic sealer from the gel.
 - (2) Assemble the gel and pre-run at 60V for about 30 min using pre-chilled 0.5X TBE as the running buffer.
 - (3) While pre-running the gel, prepare RNA samples by mixing 3 μ l of RNA loading buffer with 7 μ l (5 μ g) of total RNA, heating at 70°C for 5 min and chill on ice.
 - (4) Rinse individual wells by pipetting the buffer up and down before loading RNA samples.
 - (5) Carefully load 10 μ l RNA sample onto one well of 15% pre-run urea-polyacrylamide gel. Load 5 μ l molecular standards next to the RNA sample.
- Note: Different RNA samples can be loaded onto the gel for the detection with a single miRNA probe or a same RNA sample can be loaded onto different wells for the hybridization with different miRNA probes (see Table 1 for recommended arrangement). An empty well between two different hybridization groups is recommended for easy cutting after RNA transfer.
- (6) Run at 60V until bromophenol blue reaches approximately 3 cm away from the bottom of the gel.

2. Transfer

- (1) Disassemble the gel cast and remove one of the plates from the gel.
- (2) Transfer the gel to a glass tray filled with 0.5X TBE buffer.
- (3) Soak the membrane and filter paper in 0.5X TBE.
- (4) Assemble the transfer unit in the following order on the black side of cassette: one fiber pad, one piece of filter paper, gel, membrane and one piece of filter paper, one

(8) Dry at 42°C for 15 min.

