

HUMAN IP-10 ELISA

Product Data Sheet

Cat. No.: RGP019R

For Research Use Only

CONTENTS

1.	INTENDED USE	3
2.	INTRODUCTION	3
3.	PRINCIPLE OF THE METHOD	3
4.	REAGENTS PROVIDED AND RECONSTITUTION	4
5.	MATERIAL REQUIRED BUT NOT PROVIDED	4
6.	SAFETY	5
7.	PROCEDURAL NOTES/LAB. QUALITY CONTROL	5
8.	SPECIMEN COLLECTION, PROCESSING AND STORAGE	6
9.	PREPARATION OF REAGENTS	6
10.	ASSAY METHOD	7
11.	SUGGESTED PLATE SCHEME	9
12.	DATA ANALYSIS	10
13.	LIMITATIONS OF THE PROCEDURE	11
14.	PERFORMANCES AND CHARACTERISTICS	11
15.	REFERENCES	12
16.	ASSAY PROCEDURE SUMMARY	13

➤➤ This kit is manufactured by:
BioVendor – Laboratorní medicína a.s.

➤➤ Use only the current version of Product Data Sheet enclosed with the kit!

1. INTENDED USE

The IP-10 ELISA is to be used for research use only quantitative determination of human interferon gamma inducible 10 kD protein, in cell culture supernatants buffered solutions, serum and plasma samples or other body fluids. **This kit has been configured for research use only and is not to be used in diagnostic procedures.**

2. INTRODUCTION

IP-10 (Interferon-gamma inducible Protein 10 kDa) also known as CXCL10, is secreted by several cell types in response to IFN γ and LPS.

These cell types include monocytes, endothelial cells and fibroblasts.(1).The gene for IP-10 is located on chromosome 4 in a cluster among several other cytokines and encodes a 98 amino acid precursor protein.(1).

IP-10 has been attributed to several roles, such as chemoattraction for monocytes and T cells (but not for neutrophils), inhibition of bone marrow colony formation and angiogenesis, promotion of T cells adhesion molecule expression (2) (3).

IP-10 shares a common receptor, CXCR3, with the chemokine MIG, but has been shown to play a distinct role in host defense in infections. (4)

IP-10 expression has been associated with HIV infection (5), is involved in inflammatory skin disease (6) and other allergic diseases; it appears in inflammation of the nervous system and in Alzheimer's disease (astrocytes expressing IP-10 are commonly associated with senile plaques) (7).

3. PRINCIPLE OF THE METHOD

The IP-10 Kit is a solid phase sandwich Enzyme Linked-Immuno-Sorbent Assay (ELISA). A monoclonal antibody specific for IP-10 has been coated onto the wells of the microtiter strips provided. Samples, including standards of known IP-10 concentrations and unknowns are pipetted into these wells.

During the first incubation, the IP-10 antigen is added to wells. After washing, a biotinylated monoclonal antibody specific for IP-10 is incubated. Then the enzyme (streptavidin-horse radish peroxidase) is added. After incubation and washing to remove all unbound enzyme, a substrate solution of the bound enzyme is added to induce a coloured reaction product. The intensity of this coloured product is directly proportional to the concentration of IP-10 present in the samples.

4. REAGENTS PROVIDED AND RECONSTITUTION

REAGENTS (Store at 2-8°C)	COLOUR CODE	Quantity	State
Antibody Coated Microtiter Strips		96 wells	Ready to use
Plastic cover		2	
Standard: 200 pg/ml	Yellow	2 vials	Reconstitute with the volume of standard diluent indicated in the Quality Control Sheet. (See Reagents Preparation on page 2)
Standard Diluent buffer	Black	1 vial (25 ml)	10X concentrated. Dilute in distilled water
Biotinylated anti-IP-10	Red	1 vial (0.4 ml)	Dilute in biotinylated antibody diluent
Biotinylated Antibody Diluent	Red	1 vial (7 ml)	Ready to use
Streptavidin-HRP		2 vials (5 µl)	0.5 ml of HRP-Diluent before further dilutions
HRP Diluent	Red	1 vial (23 ml)	Ready to use
Washing Buffer	White	1 vial (10 ml)	200X concentrated. Dilute in distilled Water
Chromogen TMB :		1 vial (11 ml)	Ready to use
H ₂ SO ₄ : Stop Reagent	Black	1 vial (11 ml)	Ready to use

5. MATERIAL REQUIRED BUT NOT PROVIDED

- Distilled water.
- Pipettes : 10 µl, 50 µl, 100 µl, 200 µl and 1000 µl.
- Vortex mixer and magnetic stirrer.

6. SAFETY

- For research use only.
- Avoid any skin contact with H_2SO_4 and TMB. In case of contact, wash thoroughly with water.
- Do not eat, drink, smoke or apply cosmetics where kit reagents are used.
- Do not pipette by mouth.

7. PROCEDURAL NOTES/LAB. QUALITY CONTROL

1. When not in use, kit components should be stored refrigerated or frozen as indicated on vials or bottles labels. All reagents should be warmed to room temperature before use. Lyophilised standards should be discarded after use.
2. Once the desired number of strips has been removed, immediately reseal the bag to protect the remaining strips from deterioration.
3. Cover or cap all reagents when not in use.
4. Do not mix or interchange reagents between different lots.
5. Do not use reagents beyond the expiration date of the kit.
6. Use a clean disposable plastic pipette tip for each reagent, standard, or specimen addition in order to avoid cross-contamination; for the dispensing of H_2SO_4 and substrate solution, avoid pipettes with metal parts.
7. Use a clean plastic container to prepare the washing solution.
8. Thoroughly mix the reagents and samples before use by agitation or swirling.
9. All residual washing liquid must be drained from the wells by efficient aspiration or by decantation followed by tapping the plate forcefully on absorbent paper. Never insert absorbent paper directly into the wells.
10. The TMB solution is light sensitive. Avoid prolonged exposure to light. Also, avoid contact of the TMB solution with metal to prevent colour development. Warning TMB is toxic avoid direct contact with hands. Dispose off properly.
11. If a dark blue colour develops within a few minutes, this indicates that the TMB solution has been contaminated and must be discarded. Read absorbances rapidly after completion of the assay.
12. When pipetting reagents, maintain a consistent order of addition from well-to-well. This will ensure equal incubation times for all wells.
13. Respect incubation times described in the assay procedure.

8. SPECIMEN COLLECTION, PROCESSING AND STORAGE

Cell culture supernatants - Remove particulates and aggregates by spinning at approximately 1000 g for 10 min.

Cell lysats - After spinning at approximately 400 g for 5 min, remove the supernatant and wash once again with PBS. Suspend cells in a cold lysis buffer. After 30 min of incubation, carefully remove the supernatant after spin at 10000 g for 10 min at 4°C. Store at -70°C.

Serum - Avoid any unintentional stimulation of the cells by the procedure. Use pyrogen/endotoxin free collecting tubes. Serum should be removed rapidly and carefully from the red cells after clotting. For that, after clotting, centrifuge at approximately 1000 g for 10 min and remove serum.

Storage - If not analysed shortly after collection, samples should be aliquoted (250-500 µl) to avoid freeze-thaw cycles and stored frozen at -70°C. Avoid multiple freeze-thaw cycles of frozen specimens. When possible, avoid use of badly hemolyzed or lipemic sera. If large amounts of particles are present, this should be removed prior to assay by centrifugation or filtration.

Recommendation: Do not thaw by heating at 37°C or 56°C. Thaw at room temperature and make sure that sample is completely thawed and homogeneous before assaying.

9. PREPARATION OF REAGENTS

9.1 Standard diluent Buffer 10X concentrated

Dilute 10 times in distilled water.

9.2 Standards

Standard have to be reconstituted with the volume of standard buffer diluent indicated in the Quality Control Sheet. This reconstitution produces a stock solution of 200 pg/ml IP-10. Serial dilutions of standard must be made before each assay and cannot be stored.

9.3 Samples: serum or plasmas:

It is recommended to dilute 1:2 (see Assay method e))

9.4 Dilution of biotinylated anti - IP-10

Preparation immediately before use is recommended. Dilute the biotinylated anti-IP-10 with the biotinylated antibody diluent in a clean glass vial according to the number of wells to be used. See the next table for volumes to pipette. Extemporaneous preparations are recommended.

Number of Wells used	Biotinylated Antibody (µl)	Biotinylated Antibody Diluent (µl)
16	40	1060
24	60	1590
32	80	2120
48	120	3180
96	240	6360

9.5 Dilution of Streptavidin-HRP

Add 0.5 ml of HRP diluent to a 5 µl vial of Streptavidin-HRP. DO NOT KEEP THIS DILUTION FOR FURTHER EXPERIMENTS.

Dilute immediately before use. Following the number of wells to be used, further dilutions of Streptavidin-HRP should be made with HRP diluent in a clean glass vial : see hereafter the table for volumes to pipette.

Number of Wells	Streptavidin-HRP(µl)	Strep-HRP Diluent (ml)
16	30	2
24	45	3
32	60	4
48	75	5
96	150	10

9.6 Washing Buffer 200X concentrate

Dilute 200 times in distilled water.

10. ASSAY METHOD

- a) Before use, mix all reagents thoroughly without making foam.
- b) Determine the number of microwell strips required to test the desired number of samples, plus appropriate number of wells needed for running blanks and standards. Each sample, standard and blank should be assayed **in duplicate** . Remove sufficient microwell strips from the pouch.
- c) Add 100 µl of standard diluent (see preparation of reagents) to standard wells B1, B2, C1, C2, D1, D2, E1, E2, F1, F2. Reconstitute standard vial with the appropriate volume as described in the chapter reagents preparation. Pipet 200 µl of standard into wells A1 and A2 (see Plate Scheme below). Transfer 100 µl from A1 and A2 to B1 and B2 wells. Mix the contents by repeated aspirations and ejections. Take care not to scratch the inner surface of microwells. Repeat this procedure from the wells B1, B2 to wells C1, C2 and from wells C1, C2 to D1, D2 and so on creating two parallel rows of IP-10 standard dilutions ranging from 6.25 to 200 pg/ml. Discard 100 µl from the content of the last microwells used (F1, F2). Alternatively these dilutions can be done in separate tube and diluted standard pipetted directly into wells.

- d) Add 100 µl of standard diluent to the blank wells (G1-G2).
- e) Add 50 µl of sample and 50 µl of standard diluent to sample wells.
- f) Cover with a plate cover and incubate for 2 hours at room temperature (18°C - 25°C).
- g) Remove the cover and wash the plate as follows:
 - 1) aspirate the liquid from each well;
 - 2) dispense 0.3 ml of washing solution into each well;
 - 3) aspirate again the content of each well;
 - 4) Repeat steps 2) and 3) two times.
- h) Preparation of biotinylated anti-IP-10: (see preparation of reagents).
- i) Add 50 µl of diluted biotinylated anti-IP-10 to all wells.
- j) Cover and incubate 1 hour at room temperature.
- k) Wash as described in point g)
- l) Prepare HRP solution just before use: (see preparation of reagents).
- m) Dispense 100 µl of HRP solution into all wells, including the blank wells. Put back the cover.
- n) Incubate the microwell strips at room temperature for 30 minutes.
- o) Remove plate cover and empty wells. Wash microwell strips according to point g). Proceed immediately to the next step.
- p) Pipette 100 µl of ready-to-use TMB substrate solution into all wells, including the blank wells and incubate in the dark for 5-15 min minutes at room temperature. Avoid direct exposure to light by wrapping the plate in aluminium foil. Incubation time of the substrate solution is usually determined by the ELISA reader performances: many ELISA readers record absorbance only up to 2.0 O.D. The O.D. values of the plate should be monitored and the substrate reaction stopped before positive wells are no longer properly readable.
- q) The enzyme-substrate reaction is stopped by quickly pipetting 100 µl of H₂SO₄ : stop reagent into each well, including the blank wells, to completely and uniformly inactivate the enzyme. Results must be read immediately after the addition of H₂SO₄ : stop reagent.
- r) Read absorbance of each well on a spectrophotometer using 450 nm as the primary wavelength and optionally 620 nm (610 nm to 650 nm is acceptable) as the reference wavelength.


11. SUGGESTED PLATE SCHEME

Standard Concentrations pg/mL		Sample wells											
1	2	3	4	5	6	7	8	9	10	11	12		
A	200												
B	100												
C	50												
D	25												
E	12,5												
F	6,25												
G	Blank												
H													

12. DATA ANALYSIS

Generate a linear standard curve by plotting the average absorbance on the vertical axis versus the corresponding IP-10 standard concentration on the horizontal axis. The amount of IP-10 in each sample is determined by extrapolating OD values to IP-10 concentrations using the standard curve.

Multiply results by the appropriate dilution factor : if samples have been diluted according the recommended protocol, the concentration read from the standard curve must be multiplied by the dilution factor (x2)



13. LIMITATIONS OF THE PROCEDURE

Do not extrapolate the standard curve beyond the 200 pg/ml standard curve point. The dose-response is non-linear in this region and good accuracy is difficult to obtain. Concentrated samples (> 200 pg/ml) have to be tested with some dilutions with standard diluent or with your own sample buffer.

The influence of various drugs, aberrant sera (hemolyzed, hyperlipidemic, jaundiced...) has not been investigated. The rate of degradation of native IP-10 in various matrices has not been investigated.

14. PERFORMANCES AND CHARACTERISTICS

14.1 Sensitivity

The minimum detectable dose of IP-10 is 5.7 pg/ml.

This has been determined by adding 3 standard deviations to the mean optical density obtained when the zero standard was assayed 32 times in duplicates.

14.2 Precision

Intra-Assay					Inter-Assay				
Sample	n	Mean (pg/mL)	SD	CV%	Sample	n	Mean (pg/mL)	SD	CV%
A	6	60,1	0,9	1,6	A	18	56,9	2,7	4,7
B	6	32,4	0,5	1,5	B	18	29,9	2,0	6,5
C	6	19,6	0,3	1,3	C	18	17,7	1,5	8,5

14.3 Expected values

21 sera from apparently healthy donors were evaluated for the presence of IP-10 in this assay. The detected human IP-10 levels ranged between 61 and 177pg/ml with a mean level of 106.5pg/ml and a standard deviation of ± 28.6 pg/ml. The normal levels measured may vary with sample collective used.

14.4 Recovery

We obtained a recovery of 118.7% and 122.4% for IP-10 concentration when we add IP-10 in 2 different sera.

14.5 Sample stability

Aliquots of spiked serum samples were stored at -20°C and thawed up to five times : there was no significant loss of IP-10 immunoreactivity.


Aliquots of spiked serum samples were stored at -20°C, +4°C, Room temperature and 37°C and the IP-10 level was determined after 24h. There was no significant loss of IP-10 during storage under above conditions.

15. REFERENCES


1. Luster, A.D. et al. (1985) Nature 315:672.
2. Dufour et al. (2002) J.Immunol.168
3. Angiolillo et al. (1995) J.Exp.Med.182 :155.
4. Both et al. (2002) Biochemistry .41
5. Kolb et al. (1999) J.Neuroimmunol.108
6. Sebastiani et al. (2002) Arch.Dermatol.Res. 293
7. Strieter et al. (1995) Shock. 155-60


16. ASSAY PROCEDURE SUMMARY

Total procedure length 3 h 45 mn



NOTES





Gentaur Molecular Products
Voortstraat 49
1910 Kampenhout, Belgium
<http://www.gentaur-worldwide.com>