



# Human IGF-1 ELISA

**(Insulin Like Growth Factor-1)**


**Cat. No.: RMEE20**



## INTRODUCTION

Insulin-like growth factors (IGF) -1 and -2 play a pivotal role in regulating the proliferation, differentiation and specific functions of many cell types (1-3). IGF-1 is identical with Somatomedin C (Sm-C) (4) and has a molecular weight of 7649 daltons (5). Its major regulators are growth hormone (GH) and nutrition (6), although its production in specific tissues is affected by a multitude of tropic hormones and other peptide growth factors. In contrast to many other peptide hormones, IGFs are avidly bound to specific binding proteins (IGFBP). The seven classes of IGFBPs which are known at present (7,8,22) either bind IGF-1 and IGF-2 with similar affinities or show a preference for IGF-2 (9,10).

A major problem of IGF-1 measurement results from the interference of IGFBPs in the assay. Direct determinations in untreated serum samples (11) give false values because of the extremely slow dissociation of the IGF-1/IGFBP-3 complexes during the assay incubation. Depending on the ratio IGF-1 to IGFBP the following errors may occur (see also Figure 1):



**Figure 1.** Interference of IGFBP in IGF-1 measurements. Known concentrations of IGF-1 were assayed in the presence of 0.5 ng (left) or 5 ng (right) hIGFBP-3 by a conventional (□) and by an IGFBP-blocked assay (\*).

Therefore, various techniques were applied to physically separate IGF-1 from its binding proteins before measurement, including (a) size exclusion chromatography under acidic conditions, (b) solid-phase extraction and (c) acid-ethanol extraction (2,12,13). These techniques, however, are either inconvenient or time-consuming or give incomplete and not-reproducible recoveries. The most widely used method is the acid-ethanol extraction (13,14) with a recovery of only 70-80 % of IGFBP-bound IGF-1 as a result of co-precipitation. The absolute results of such an extraction are therefore false low (15). The extraction removes the IGFBPs only insufficiently and leads to reduction in sensitivity of the assay due to pre-dilution of the samples by the extraction procedure. Furthermore, the remaining IGFBP may still interfere in the assay. In addition, the acid-ethanol extraction is ineffective in specimens other than serum or plasma (e.g. cell culture media), in which determination of IGF-1 is already difficult enough due to the fact that IGFBPs are frequently present at large excess.

To avoid these difficulties, an uncomplicated assay was developed, in which special sample preparation is not required before measurement.

## CLINICAL SIGNIFICANCE


There are apart from GH, a number of variables that influence serum IGF-1. Decreased levels are found in states of malnutrition/malabsorption, hypothyroidism, liver disease, untreated diabetes mellitus, chronic inflammatory disease (1,6), malignant disease or polytrauma. High levels, on the other hand, are likely to be present in precocious puberty or obesity. Crucially important to the correct interpretation of IGF-1 measurements is the relationship between age and IGF-1 levels. It is certainly inadequate to use a common cut-off point to define "normal" levels for all age groups, particularly in children and adolescents.

Due to its GH-dependence, determination of serum IGF-1 was shown to be a useful tool in diagnosis of growth disorders, especially with regard to GH deficiency (GHD) or acromegaly (6,16-19,23,24). The major advantage of IGF-1 determination compared to GH determination is its stable circadian concentration; therefore a single measurement is sufficient. Hence IGF-1 determination should be the first in a series of laboratory test. Clearly normal levels would then rule out disturbances of the GH-IGF-1-axis. Low levels, i.e. close to or below the age-related 5th percentile would indicate the necessity of further diagnostic efforts. Subnormal levels of IGF-1 would be evidence for reduced GH secretion, if other causes of low serum IGF-1 (e.g. malnutrition or impaired liver function) can be ruled out. For differentiation of healthy short children without GH deficiency and children with "classical" GH deficiency, the 0.1st percentile proved to be an appropriate cut-off point, especially after the age of eight. However, IGF-1 levels of short children not suffering from GHD may nevertheless lay between the 0.1st and 5th percentile (19). In contrast, acromegaly is characterized by pathologically elevated IGF-1 levels, which apparently reflect the severity of the disease better than GH-levels (17,18,20).

## INTENDED USE

An enzyme immunoassay for the quantitative *in vitro diagnostic* measurement of levels of IGF-1 in serum or plasma. Human IGF-1 measurements are used in the diagnosis and treatment of growth disorders involving the anterior lobe of the pituitary gland.

## PRINCIPLE



**Figure 2.:** Principle of the IGFBP-blocked IGF-1 ELISA

In order to dissociate IGF-1 from the IGFBPs, the samples must be diluted in an acidic buffer (Figure 2). The diluted samples are then pipetted into the assay wells. The IGF-1 antiserum is dissolved in a buffer, which is able to neutralize the acidic samples. After the IGF-1 antibody solution has neutralized the samples, the present excess IGF-2 occupies the IGF-binding sites of the binding proteins, thus allowing the measurement of the resulting free IGF-1. With this method, the IGFBPs are not removed, but their function and therefore their interference in the assay is neutralized. Due to the extremely low cross-reactivity of the IGF-1 antibody with IGF-2, the excess of IGF-2 does not disturb the interaction of the first antibody with IGF-1.

The ELISA for IGF-1 RMEE20 is a so-called Sandwich-Assay using two specific and high affinity antibodies. The IGF-1 in the samples binds to the first antibody coated on the microtiterplate, the second specific anti-IGF-1-antibody binds

in turn to the immobilised IGF-1. The second antibody is biotinylated, the subsequently incubated Streptavidin-Peroxidase-Enzyme Conjugate will bind to it, and thus in the final substrate incubation step colour development will be catalysed quantitatively depending on the IGF-1-level of the samples. The Standards of the ELISA RMEE20 are prepared from recombinant IGF-1 in concentrations of 2, 5, 15, 30 and 50 ng/ml.

## PERFORMANCE CHARACTERISTICS AND VALIDATION

### Sensitivity

The analytical sensitivity of the enzyme immunoassay for IGF-1 yields **0.09 ng/mL** 2x SD of zero standards in 17-fold determination.

### Specificity

The following materials have been evaluated for cross reactivity. 200 ng/ml solutions of each substance have been analysed in this Enzyme Immunoassay.

	IGF-2	Insulin	Proinsulin	C-Peptide
Reactivity [%]	0.055	0.045	0.025	0.02

## Reproducibility

### Intra-Assay-Variation

	Number of determinations	Mean value (ng/ml)	Standard deviation (ng/ml)	CV (%)
Sample 1	18	144.8	9.63	6.65
Sample 2	18	140.79	7.15	5.08
Sample 3	18	138.02	7.86	5.69

### Inter-Assay-Variation

	Mean value (ng/ml)	Standard deviation (ng/ml)	CV (%)
Sample 1	174	11.79	6.79
Sample 2	494	11.11	2.25
Sample 3	142	8.68	6.11

## Linearity

Dilution:	Sample 1 (calculated, ng/ml)	Dilution:	Sample 2 (calculated, ng/ml)
1:10	137.2	1:10	439.1
1:20	133.5	1:20	500.2
1:40	133.6	1:40	499.2
1:80	134.6	1:80	490.5
1:160	134.4	1:160	494.5
1:320	135.7	1:320	526.4
		1:640	463.7
<b>AV / 1SD / VC%</b>	134.8/ 1.4/ 1.04	<b>AV / 1SD / VC%</b>	487.6/ 28.2/ 5.79

AV = Average Value , SD = Standard Deviation

## Recovery

The **recovery** of the recombinant hIGF-1 yielded in a buffer matrix 100%. In different human-sera the recovery was on average 98.67% (n=5) of the hypothetical expected amount.

#### Exemplary data

	Serum 1	Serum 2	Serum 3
IGF-1 added	200 ng/ml	400 ng/ml	400 ng/ml
Recovery [%]	103.4	101.1	101.9

#### SPECIMEN COLLECTION, PREPARATION, AND STORAGE

The stability of IGFBP-bound IGF-1 makes sample preparation simple. Blood samples may be taken at any time of the day. Whole blood should be processed within a few hours and stored frozen at  $-4^{\circ}\text{F}$  ( $-20^{\circ}\text{C}$ ) until measurement. IGF-1 levels are usually not affected by improper handling or storage. They remain stable over several days in normal and in various clinical situations even under conditions of high temperature  $98.6^{\circ}\text{F}$  ( $37^{\circ}\text{C}$ ). Avoid repeated freezing and thawing cycles, although IGF-levels in normal sera remained unchanged after 10 cycles. Frozen samples are stable over many years. Samples may also be freeze-dried without suffering any loss of activity.

#### Sample Preparation

Serum samples as well as Heparin-, EDTA- and Citrat-Plasma samples are suited. Possible dilution of the sample by the anticoagulant must be considered.

Serum or plasma samples should be diluted depending on the expected values 1:10 -1:50 with **Sample Buffer PP**. Generally, in case of serum or plasma specimens a dilution of **1:21** is very well suited. IGF-1 concentration in other body fluids or cell culture supernatants could differ strongly.

##### Suggestion for dilution protocol:

Pipette 200  $\mu\text{l}$  **Sample Buffer PP** in PE/PP-Tubes (application of a multi-stepper is recommended in larger series); add 10  $\mu\text{l}$  serum or plasma samples (dilution 1:21). After mixing use  $2 \times 20 \mu\text{l}$  of this dilution in the assay.

## REAGENTS PROVIDED

1)	MTP	<b>Microtiter plate</b> , ready for use, with 96 wells, divided up in 12 strips with 8 wells separately breakable, coated with anti-human IGF-1 antibody, packed in a laminate bag.
2)	CAL	<b>Standards A-E</b> , lyophilised, contain recombinant human IGF-1. Standard values are <b>between 2-50 ng/ml</b> (2, 5, 15, 30 and 50 ng/ml) IGF-1 and have to be reconstituted in <b>500 µl (each) in Sample Buffer PP</b> . After using store the reconstituted standards in the original flasks as soon as possible at $-20^{\circ}\text{C}$ ( $-4^{\circ}\text{F}$ ). When using the standards anew, please thaw them rapidly but gently (no temperature rise over the room temperature and no powerful vortexing), 3 of these freezing-thawing cycles showed no influence on the assay. 20 µl per well are used in the assay.
3)	BUF PP	<b>Sample Buffer PP, 25 ml</b> , ready for use, please use for reconstitution of Standards and Controls and for dilution of Samples and Controls
4)	Control	<b>Control Serum KS1 and KS2, 500 µl</b> , lyophilised, contain human serum and has to be reconstituted in <b>500 µl Sample Buffer PP</b> . The reconstituted Control Sera must be stored in the original flask as soon as possible at $-20^{\circ}\text{C}$ ( $-4^{\circ}\text{F}$ ) after using. When using anew, please thaw them rapidly but gently (no temperature rise over the room temperature and no powerful vortexing), 3 of these freezing-thawing cycles showed no influence on the assay. The IGF-1 target value concentrations and the respective ranges are given on the vial labels. The dilution of the Control Sera should be according to the dilution of the respective samples.
5)	Ab	<b>Antibody Conjugate AK, 9 ml</b> , ready for use, contains the biotinylated anti-IGF-1 antibody. Use 80 µl for each well in the assay. <b>ATTENTION: READY FOR USE!</b>
6)	CONJ	<b>Enzyme Conjugate EK, 12 ml</b> , ready for use, contains horseradish-peroxidase conjugate to streptavidin, use 100 µl for each well in the assay. <b>ATTENTION: READY FOR USE!</b>
7)	WASHBUF 20x	<b>Washing Buffer (WP), 50 ml, 20 X concentrated solution. Washing Buffer (WP), 50 ml, 20X concentrated solution.</b> Washing Buffer (WP) has to be diluted 1:20 with distilled or demineralised water before use (e.g. add the complete contents of the flask (50 ml) into a graduated flask and fill up with A.dest. to 1000 ml). Attention: After dilution the Washing Buffer is only 4 weeks stable, dilute only according to requirements.
8)	SUBST	<b>Substrate (S), 12 ml</b> , ready for use, horseradish-peroxidase-(HRP)-substrate, stabilised $\text{H}_2\text{O}_2$ Tetramethylbencidine.
9)	$\text{H}_2\text{SO}_4$	<b>Stopping Solution (SL), 12 ml</b> , ready for use, 0.2 M sulphuric acid, Caution acid!
10)		Sealing tape for covering of the microtiter plate, 2 x, adhesive.

## MATERIALS REQUIRED BUT NOT PROVIDED

- Distilled or demineralised water for dilution of the Washing Buffer WP
- Micropipettes and multichannel pipettes with disposable plastic tips
- Vortex-mixer
- Microtiter plate Shaker (350 rpm)
- Microtiter plate washer (recommended)
- Microtiter plate reader ("ELISA-Reader") with filter for 450 and  $\geq 590$  nm.
- Polyethylen PE/Polypropylen PP tubes for dilution of samples

## REAGENT PREPARATION

**Bring all reagents to room temperature 68-77°F (20 - 25°C) before use.**

Wash Solution: Add distilled/deionized water to the complete contents of the (20x) concentrated Wash Solution (50 ml) into a graduated flask and fill to a final volume of 1000 ml, or dilute 1:20 only according to requirements. The diluted Washing Buffer is stable for max. 4 weeks at 35.6-46.4°F (2-8°C).

For the reconstitution of the lyophilised components (Standards A - E and Control Sera KS1 & KS2) Sample Buffer PP should be used. It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer.

The incubation steps should be performed at mean rotation frequency of a particularly suitable microtiter plate shaker. We are recommending 350 rpm. Due to certain technical differences deviations may occur, in case the rotation frequency must become adjusted. Insufficient shaking may lead to inadequate mixing of the solutions and thereby to low optical densities, high variations and/or false values, excessive shaking may result in high optical densities and/or false values.

Proper washing is of basic importance for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided washing buffer diluted to usage concentration. Washing volume per washing cycle and well must be 300 µl at least.

The danger of handling with potentially infectious material must be taken into account.

When using an automatic microtiter plate washer, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non fuzzy absorbent tissue.

Manual washing is an adequate alternative option. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamically swinging out the microtiter plate over a basin. If aspirating devices are used, care has to be taken that the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non fuzzy absorbent tissue.



## STORAGE CONDITIONS

The **microtiterplate wells** and all undiluted reagents are stable until the expiry date, if stored in the dark at 35.6-46.4°F (2-8°C).

Store the microtiterplate wells sealed together with the desiccant at 35.6-46.4°F (2-8°C).

The **Substrate Solution S**, stabilised H<sub>2</sub>O<sub>2</sub>-Tetramethylbenzidine, is photosensitive – store and incubate in the dark.

Reconstituted components (**Standards A – E** and **Control Sera KS1 & KS2**) should be stored at –4°F (-20°C). Freezing extends the expiry at least 2 months. When using the Standards or Control Sera KS1 or KS2 anew, please thaw them rapidly but gently (no temperature rise over the room temperature and no powerful vortexing), 3 of these freezing-thawing cycles showed no influence on the assay.

Precipitates, possible in buffers, should be dissolved before use through mixing and warming.

**Room temperature incubation means: Incubation at 68-77°F (20-25°C).**

## WARNINGS AND PRECAUTIONS

*For In Vitro Diagnostic Use.*

For professional use only.

**Before starting the assay, read the instructions completely and carefully.** Use the valid version of the package insert provided with the kit. Be sure that everything is understood.

Before use, all kit components should be brought to room temperature at 68-77°F (20 - 25°C). Precipitates in buffers should be dissolved before use by thorough mixing and warming. Temperature WILL affect the absorbance readings of the assay. However, correctness of the results will not be affected.

Do not mix reagents of different lots. Do not use expired reagents.

The microtiterplate contains break apart strips. Unused wells must be stored at 35.6-46.4°F (2 - 8°C) (in the sealed foil pouch and used in the frame provided).

**Caution:** This kit contains material of human and/or animal origin. Source human sera for the Control Sera provided in this kit was tested by FDA recommended methods and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV) antibodies. No known test methods can offer total assurance of the absence of infectious agents; therefore all components and patient's specimens should be treated as potentially infectious.

**Stop solution contains 0.2 M Sulfuric Acid (H<sub>2</sub>SO<sub>4</sub>)**

R36/38 Irritating to eyes and skin

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice

S28.1 After contact with skin, wash immediately with plenty of water

S36/37 Wear suitable protective clothing and gloves.

Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.

Use separate pipette tips for each sample, control and reagent to avoid cross contamination.

Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored.

Do not pour reagents back into vials as reagent contamination may occur.

Mix the contents of the microtiterplate wells thoroughly to ensure good test results. Do not reuse microtiterplate wells.

Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.

Some reagents contain <0.01% 2-Methyl-Isothiazolin-3-one and <0.01% 5-Chloro-2-Methyl-4-Isothiazolin-3-one and 2-Methyl-Isothiazolin-3-one, as preservatives. In case of contact with eyes or skin, flush immediately with water.

**TMB-Substrate (S) contains 3,3',5,5' Tetramethylbenzidine.**

R20/21/R22 Harmful by inhalation, in contact with skin and if swallowed

R36/37/38 Irritating to eyes, respiratory system and skin

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice

S28.1 After contact with skin, wash immediately with plenty of water

S36/37 Wear suitable protective clothing and gloves

General first aid procedures:

*Skin contact:* Wash affected area thoroughly with water. Discard contaminated cloths and shoes.

*Eye contact:* In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.

*Ingestion:* If swallowed, wash out mouth thoroughly with water. Immediately see a physician.

Do not eat, drink or smoke in these areas.

Never pipette the materials with the mouth.

Spilled material must be wiped off immediately and should become disinfected. Clean contaminated areas and equipment with a suitable detergent.

## ASSAY PROCEDURE

**NOTES:** All determinations (**Standards**, **Control Sera** and **samples**) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

When performing the assay, the **Standards**, **Control Sera** and the samples should be pipette as fast as possible (e.g., <15 minutes). To avoid distortions due to differences in incubation times, the **Enzyme Conjugate** as well as the succeeding **Substrate Solution** should be added to the plate in the same order and in the same time interval as the samples. **Stop Solution** should be added to the plate in the same order as the **Substrate Solution**.

1. Add **80 µl Antibody Conjugate AK** in **all** wells used
2. Pipette in positions A1/2 **20 µl Sample Buffer PP**
3. Pipette in positions B1/2 **20 µl of the Standard A (2 ng/ml)**  
Pipette in positions C1/2 **20 µl of the Standard B (5 ng/ml)**,  
Pipette in positions D1/2 **20 µl of the Standard C (15 ng/ml)**,  
Pipette in positions E1/2 **20 µl of the Standard D (30 ng/ml)**,  
Pipette in positions F1/2 **20 µl of the Standard E (50 ng/ml)**.

To control the correct accomplishment of the assay **20 µl** of the 1:21(or in respective dilution ratio of the samples) in **Sample Buffer PP** diluted **Control Sera KS1&KS2** can be pipetted in positions G1/2 and H1/2.

Pipette **20 µl** each of the diluted sample (e.g. dilute 1:21 with Sample buffer **PP**) in the rest of wells, according to your requirements.

4. Cover the wells with sealing tape and incubate the plate for **1 hour** at **room temperature** 68-77°F (20-25°C) (shake at 350 rpm).
5. After incubation aspirate the contents of the wells and wash the wells 5 times **300 µl Washing Buffer WP** / well.
6. Following the last washing step pipette **100 µl** of the **Enzyme Conjugate EK** in each well.
7. Cover the wells with sealing tape and incubate the plate for **30 minutes** at **room temperature** 68-77°F (20-25°C) (shake at 350 rpm).
8. After incubation wash the wells 5 times with **Washing Buffer WP** as described in step 5.

9. Pipette **100 µl** of the **Substrate Solution S**.
10. Incubate the plate for **15 minutes in the dark** at **room temperature** 68-77°F (20-25°C).
11. Stop the reaction by adding **100 µl Stopping Solution SL** to all wells.

Measure the absorbance within 30 minutes at **450 nm (Reference filter  $\geq$  590 nm)**.

## QUALITY CONTROL

The handling of potentially infectious material must comply with Good Laboratory Practice (GLP). GLP requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance. The test results are only valid if the test has been performed following the instructions. Moreover the user must strictly adhere to the rules of GLP or other applicable federal, state or local standards/laws. All standards and kit controls must be found within the acceptable ranges as stated on the QC Certificate. If the criteria are not met, the run is not valid and should be repeated. Each laboratory should use known samples as further controls.

## CALCULATION OF RESULTS

### Assay Characteristics and Validation

The standards are derived from recombinant hIGF-1 devoid of methIGF-1 or IGF-1 variants with mismatched disulfide bonds, i.e. this recombinant IGF-1 is identical to the major authentic IGF-1 form in blood. The Mediagnost IGF-1 enzymeimmunoassay is calibrated against the International Reference Standard preparation of IGF-1, **WHO NIBSC Code 02/254** (25, 26)

For the evaluation of the assay it is required that the absorbance values of the blank should be below 0.25 and the absorbance of standard E should be above 1.00. Samples, which yield higher absorbance values than **Standard E**, are beyond the standard curve. For reliable determinations these samples should be re-tested with a higher dilution.

The standards provided contain the following concentration of recombinant hIGF-1:

Standard	A	B	C	D	E
ng/ml	2	5	15	30	50

- 1) Calculate the mean absorbance (MA) value for the blank from the duplicated determination (well A1/A2).
- 2) Subtract the means absorbance of the blank from the mean absorbances of the standards and of the samples.
- 3) Plot the standard concentrations on the x-axis versus the mean value of the absorbance of the standards on the y-axis.

- 4) Recommendation: Calculation of standard curve should be done by using a computer programme, because the curve is in general (without respective transformation) not ideally described by linear regression. **A higher-grade polynomial, or four parametric logistic (4-PL) curve fit or non-linear regression** are usually suitable for the evaluation (as might be spline or point-to-point alignment in individual cases).
- 5) The IGF-1 **concentration of the sample in ng/ml** can be calculated by **multiplication** with the respective **dilution factor**.

#### Example of Typical Standard Curve

The following data is for demonstration only and cannot be used in place of data generation at the time of assay.

Standard	A	B	C	D	E
ng/ml	2	5	15	30	50
extinction	0.07	0.18	0.59	1.06	1.58

## EXPECTED NORMAL VALUES

IGF-1 levels are highly age-dependent in children, less so in adults until the age of about 60. The normal ranges in various age groups, which are log-normally distributed, are given in Table 2 by percentiles. Between 8 and 19 years of age, values are given for boys and girls separately, because the pubertal peak usually occurs approximately 2 years earlier in girls. A graphic presentation is shown in Figures 3, 4 and 5. A major problem for the interpretation of IGF-1 values arises from the fact that short stature is often due to developmental delay rather than any metabolic or endocrine disorder (constitutional delay of growth and adolescence). The sharp rise in IGF-1 levels during puberty may therefore cause some uncertainty as to whether or not it would be appropriate to relate measured values to chronological age. It is recommended to take the pubertal stage into account (Table 1 and Figure 6) to get a more complete picture of this situation.

**Table 1:** Normal range of serum IGF-1 levels at different pubertal stages according to Tanner. Because no significant difference between boys and girls is observed, both sexes are combined. Only children and adolescents between 7 and 17 years of age are included.

Percentile				
Pubertal Stage	0.1th	5th	50th	95th
1	61	105	186	330
2	85	156	298	568
3	113	196	352	631
4	171	268	431	693
5	165	263	431	706

**Table 2:** Serum levels of IGF-1 in healthy subjects at various ages. Individuals between 8 and 19 years of age were classified according to gender, as the pubertal peak occurs almost 2 years earlier in girls than in boys.

<b>Percentile</b>														
Age	0.1	1	5	10	20	30	40	50	60	70	80	90	95	99
0-2 y.	13	20	28	34	43	50	58	66	75	87	102	128	156	220
2-4 y.	20	29	40	48	59	68	77	87	98	111	129	159	189	260
4-6 y.	26	36	50	59	73	85	96	108	122	138	160	196	233	320
6-7 y.	34	46	62	72	87	99	111	124	138	155	176	212	248	332
7-8 y.	45	60	78	90	107	121	134	148	163	181	205	243	281	364
8-9 y.	54	71	90	102	119	133	146	160	175	192	214	250	284	362
		75	99	115	137	156	174	193	214	239	271	324	376	496
9-10 y.	63	82	102	115	133	148	162	176	191	209	232	269	304	379
		89	114	130	152	170	187	205	224	247	276	323	369	469
10-11 y.	77	96	117	130	148	162	176	189	203	220	241	274	305	370
		106	134	153	178	199	219	239	261	287	321	374	426	539
11-12 y.	85	106	129	144	163	179	194	209	225	244	267	304	339	413
		123	160	185	220	248	276	305	337	374	424	503	581	758
12-13 y.	88	112	141	159	184	204	223	243	264	289	321	371	419	525
		155	201	231	274	309	342	377	415	460	519	614	707	914
13-14 y.	111	143	179	203	235	261	286	311	339	371	412	477	540	677
		207	256	287	329	364	395	428	463	504	556	637	716	884
14-15 y.	140	182	229	260	303	337	370	404	441	484	539	625	691	896
		236	284	314	353	385	414	443	474	510	556	628	713	832
15-16 y.	176	221	269	299	340	372	402	433	466	504	552	626	697	849
		231	279	309	350	382	412	442	474	512	559	632	700	845
16-17 y.	178	221	267	296	335	366	395	424	455	491	537	607	673	814
		225	270	298	336	366	394	422	452	486	530	597	660	792
17-18 y.	173	207	243	265	294	317	337	358	380	405	436	484	527	618
		210	246	268	297	320	341	362	384	409	441	488	533	624
18-19 y.	167	201	235	256	285	307	327	347	368	393	423	469	512	600
		199	233	254	281	302	322	341	362	385	414	458	499	583
19-20 y.	158	189	220	240	265	285	304	322	341	363	391	433	471	550
20-30 y.	72	92	115	130	150	167	182	198	215	235	261	302	340	425
30-40 y.	68	87	109	123	142	158	173	188	204	223	248	287	324	404
40-50 y.	64	82	103	116	135	150	164	178	194	212	235	272	310	385
50-60 y.	60	77	97	110	127	142	155	169	184	201	224	260	292	369
60-70 y.	55	72	91	103	120	134	147	161	176	193	215	251	282	362
70-80 y.	25	35	47	55	67	78	88	98	110	124	142	173	207	276
>80 y.	21	30	40	47	58	67	76	85	95	108	125	153	184	245




Fig. 3: Age-dependent normal range of serum IGF-I levels in girls




Fig. 4: Age-dependent normal range of serum IGF-I levels in boys




Fig. 5: Age-dependant normal range of serum IGF-I levels in adults




Fig 6.: Serum IGF-1 levels in normal children and adolescents (7 to 17 years) according to Serum concentrations are given in ng/mL. Determined with IGFBP-blocked IGF-1 RIA without extraction step (Blum and Breier 1994).

## LIMITATIONS OF PROCEDURE

IGF-1 levels depend to a great degree on GH secretion. Diminished IGF-1 values, however, do not prove GH deficiency, because a number of other factors can influence the plasma concentration of IGF-1 and must therefore be taken into account in order to make a correct interpretation. IGF-1 levels decrease during fasting (more than 1 day), as a result of malnutrition, malabsorption, cachexia, impaired hepatic function, or in hypothyroidism and untreated diabetes mellitus. They may also be low in chronic inflammatory disease and malignancies. IGF-1 levels are high in states of accelerated sexual development. In clinical situations with hyperprolactinemia or in patients with craniopharyngioma, normal levels may be observed despite GH deficiency. In late pregnancy, IGF-1 levels are moderately elevated.

**Table 5:** Serum levels of IGF-1 in healthy subjects at various ages. Individuals between 8 and 19 years of age were classified according to gender, as the pubertal peak occurs almost 2 years earlier in girls than in boys.

	Percentiles													
Age	0,1	1	5	10	20	30	40	50	60	70	80	90	95	99
0-2 y.	13	20	28	34	43	50	58	66	75	87	102	128	156	220
2-4 y.	20	29	40	48	59	68	77	87	98	111	129	159	189	260
4-6 y.	26	36	50	59	73	85	96	108	122	138	160	196	233	320
6-7 y.	34	46	62	72	87	99	111	124	138	155	176	212	248	332
7-8 y.	45	60	78	90	107	121	134	148	163	181	205	243	281	364
8-9 y. boys	54	71	90	102	119	133	146	160	175	192	214	250	284	362
girls	55	75	99	115	137	156	174	193	214	239	271	324	376	496
9-10 y. boys	63	82	102	115	133	148	162	176	191	209	232	269	304	379
girls	68	89	114	130	152	170	187	205	224	247	276	323	369	469
10-11 y. boys	77	96	117	130	148	162	176	189	203	220	241	274	305	370
girls	81	106	134	153	178	199	219	239	261	287	321	374	426	539
11-12 y. boys	85	106	129	144	163	179	194	209	225	244	267	304	339	413
girls	91	123	160	185	220	248	276	305	337	374	424	503	581	758
12-13 y. boys	88	112	141	159	184	204	223	243	264	289	321	371	419	525
girls	116	155	201	231	274	309	342	377	415	460	519	614	707	914
13-14 y. boys	111	143	179	203	235	261	286	311	339	371	412	477	540	677
girls	163	207	256	287	329	364	395	428	463	504	556	637	716	884
14-15 y. boys	140	182	229	260	303	337	370	404	441	484	539	625	691	896
girls	193	236	284	314	353	385	414	443	474	510	556	628	713	832
15-16 y. boys	176	221	269	299	340	372	402	433	466	504	552	626	697	849
girls	187	231	279	309	350	382	412	442	474	512	559	632	700	845
16-17 y. boys	178	221	267	296	335	366	395	424	455	491	537	607	673	814
girls	183	225	270	298	336	366	394	422	452	486	530	597	660	792
17-18 y. boys	173	207	243	265	294	317	337	358	380	405	436	484	527	618
girls	176	210	246	268	297	320	341	362	384	409	441	488	533	624
18-19 y. boys	167	201	235	256	285	307	327	347	368	393	423	469	512	600
girls	167	199	233	254	281	302	322	341	362	385	414	458	499	583
19-20 y.	158	189	220	240	265	285	304	322	341	363	391	433	471	550
20-30 y.	72	92	115	130	150	167	182	198	215	235	261	302	340	425
30-40 y.	68	87	109	123	142	158	173	188	204	223	248	287	324	404
40-50 y.	64	82	103	116	135	150	164	178	194	212	235	272	310	385
50-60 y.	60	77	97	110	127	142	155	169	184	201	224	260	292	369
60-70 y.	55	72	91	103	120	134	147	161	176	193	215	251	282	362
70-80 y.	25	35	47	55	67	78	88	98	110	124	142	173	207	276
>80 y.	21	30	40	47	58	67	76	85	95	108	125	153	184	245



## REFERENCES

1. Baxter RC. 1986 The somatomedins: insulin-like growth factors. *Adv Clin Chem.*25:49-115
2. Daughaday WH, Rotwein P. 1989 Insulin-like growth factors I and II. Peptide, messenger ribonucleic acid and gene structures, serum, and tissue concentrations. *Endocr Rev.* 10:68-91
3. Spencer EM (Ed.) 1991 *Modern Concepts of Insulin-Like Growth Factors*. New York: Elsevier.
4. Klapper DG, Svoboda ME, Van Wyk JJ. 1983 Sequence analysis of somatomedin-C: confirmation of identity with insulin-like growth factor-I. *Endocrinology.* 112:2215-2217.
5. Rinderknecht E, Humbel RE. 1978 The amino acid sequence of human insulin-like growth factor I and its structural homology with proinsulin. *J Biol Chem.* 253:2769-2276.
6. Clemmons DR, Van Wyk JJ. 1984 Factors controlling blood concentration of somatomedin C. *Clin Endocrinol Metab.* 13:113-143.
7. Ballard J, Baxter R, Binoux M, et al. 1989 On the nomenclature of the IGF binding proteins. *Acta Endocrinol (Copenh).* 121:751-752.
8. Drop SLS. 1992 Report on the nomenclature of the IGF binding proteins. *J. Clin Endocrinol Metab.* 74:1215-1216.
9. Martin JL, Baxter RC. 1986 Insulin-like growth factor binding protein from human plasma. Purification and characterization. *J Biol Chem.* 261:8754-8760.
10. Binkert C, Landwehr J, Mary JL, Schwander J, Heinrich G. 1989 Cloning, sequence analysis and expression of a cDNA encoding a novel insulin-like growth factor binding protein (IGFBP-2). *EMBO J.* 8:2497-2502.
11. Furlanetto RW, Underwood LE, Van Wyk JJ, D'Ercole AJ. 1977 Estimation of somatomedin-C levels in normals and patients with pituitary disease by radioimmunoassay. *J. Clin Invest.* 60:648-657.
12. Daughaday WH, Kapadia M, Mariz I. 1987 Serum somatomedin binding proteins: physiologic significance and interference in radioligand assay. *J Lab Clin Med.* 109:355-363.
13. Breier BH, Gallaher BW, Gluckman PD. 1991 Radioimmunoassay for insulin-like growth factor-I: solutions to some potential problems and pitfalls. *J Endocrinol.* 128:347-357.
14. Daughaday WH, Mariz IK, Blethen SL. 1980 Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites: a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol-extracted serum. *J Clin Endocrinol Metab.* 51:781-788.
15. Blum WF, Gallaher B, Ranke MB. 1992 An IGFBP-blocked IGF-I RIA that measures what it pretends to measure: IGF-I. 74th Annual Meeting of the American Endocrine Society. 293.
16. Rosenfeld RG, Wilson DM, Lee PDK, Hintz RL. 1986 Insulin-like growth factors I and II in evaluation of growth retardation. *J Pediatr.* 109:428-433.
17. Clemmons DR, Van-Wyk JJ, Ridgway EC, Kliman B, Kjellberg RN, Underwood LE. 1979 Evaluation of acromegaly by radioimmunoassay of somatomedin-C. *N.Engl J Med.* 301:1138-1142

18. Zapf J, Walter H, Froesch ER. 1981 Radioimmunological determination of insulin-like growth factors I and II in normal subjects and in patients with growth disorders and extrapancreatic tumor hypoglycemia. *J Clin Invest.* 68:1321-1330.
19. Blum WF. 1992 Insulin-like growth factors and their binding proteins. In: Ranke MB, ed. *Functional Endocrinology Diagnostics in Children and Adolescence*. Mannheim: J + J Verlag; 102-117.
20. Rieu M, Girard F, Bricaire H, Binoux M. 1982 The importance of insulin-like growth factor (somatomedin) measurements in the diagnosis and surveillance of acromegaly. *J Clin Endocrinol Metab.* 55:147-153.
21. Blum WF, Ranke MB, Bierich JR. 1988 A specific radioimmunoassay for insulin-like growth factor II: the interference of IGF binding proteins can be blocked by excess IGF-I. *Acta Endocrinol (Copenh).* 118:374-380.
22. Wilson EM, Oh Y, Rosenfeld RG (1997) Generation and characterization of an IGFBP-7 antibody: Identification of 31 kD IGFBP-7 in human biological fluids and Hs578T human breast cancer conditioned media. *J Clin Endocrinol Metab* Vol 82, 4:1301-1303.
23. Ranke MB, Schweizer R, Elmlinger MW, Weber K, Binder G, Schwarze CP, Wollmann HA (2000) Significance of basal IGF-I, IGFBP-3 and IGFBP-2 measurements in the diagnostics of short stature in children. *Horm Res* 54:60-68
24. Ranke MB, Schweizer R, Elmlinger MW, Weber K, Binder G, Schwarze CP, Wollmann HA (2001) Relevance of IGF-I, IGFBP-3 and IGFBP-2 measurements during GH treatment of GH-deficient and non-GH-deficient children and adolescents. *Horm Res* 55:155-124
25. Address NIBSC: Blance Lane, South Mimms, Potters Bar, Hertford EN 6 3 QG, Great Britain
26. Burns C, Rigsby P, Moore M, Rafferty B. (2009) The first International Standard for Insulin-like Growth Factor-1 (IGF-1) for immunoassay: Preparation and calibration in an international collaborative study. *Growth Hormone & IGF Research* (2009)
27. Blum WH, Breier BH (1994) Radioimmunoassays for IGFs and IGFBPs. *Growth Regulation* 4 (Suppl. 1):11-19
28. Ranke MB, Osterziel KJ, Schweitzer R, Schuett B., Weber K., Röbbel P, Vornwald A, Blumenstock G, Elmlinger MB (2003) Reference Levels of Insulin-Like-Growth Factor I in the Serum of Healthy Adults: Comparison of Four Immunoassays. *Clin Chem Lab Med* Vol 41(10):1329-1334

## SUMMARY – IGF-1 ELISA

Reconstitution / Dilution of Reagents		
<b>Standards A-E</b>	Reconstitution in <b>Sample Buffer PP</b>	<b>500 µl</b>
<b>Control Serum KS1</b>	Reconstitution in <b>Sample Buffer PP</b>	<b>500 µl</b>
<b>Control Serum KS2</b>	Reconstitution in <b>Sample Buffer PP</b>	<b>500 µl</b>
<b>Wash Buffer WP</b>	dilute in <b>A. dest.</b> (eg. total volume of 50 ml in a graduated flask and fill up to 1000 ml)	<b>1:20</b>
<b>Sample + Control Sera KS1 and KS2: dilute 1:21 in Sample Buffer PP, mix immediately, incubate at least for 15 min, max. 2h. Use 20 µl for each well in the assay.</b>		
Before conducting the assay equilibrate all <b>reagents to room temperature.</b>		

*Assay Procedure for Double Determinations:*

Pipette	Reagent	Position
80 µl	Antibody Conjugate <b>AK</b>	in <u>all</u> wells used
20 µl	Sample Buffer <b>PP</b> (blank)	A1 and A2
20 µl	Standard <b>A</b> (2 ng/ml)	B1 and B2
20 µl	Standard <b>B</b> (5 ng/ml)	C1 and C2
20 µl	Standard <b>C</b> (15 ng/ml)	D1 and D2
20 µl	Standard <b>D</b> (30 ng/ml)	E1 and E2
20 µl	Standard <b>E</b> (50 ng/ml)	F1 and F2
20 µl	Control Serum <b>KS1</b>	G1 and G2
20 µl	Control Serum <b>KS2</b>	H1 and H2
20 µl	Samples	following wells
Cover the wells with the sealing tape.		

**Incubation: 1 h at RT, ≥ 350 rpm**

5x 300 µl	Aspirate the contents of the wells and wash <b>5x</b> with <b>300 µl Wash Buffer WP</b>	each well
100 µl	Enzyme Conjugate <b>EK</b>	each well

**Incubation: 30 min at RT, ≥350 rpm**

5x 250 µl	Aspirate the contents of the wells and wash <b>5x</b> with <b>250 µl Wash Buffer WP</b>	each well
100 µl	Substrate <b>S</b>	each well

**Incubation: 15 min in the dark RT**

100 µl	Stop Solution <b>SL</b>	each well
Measure the absorbance within <b>30 min at 450 nm</b> with <b>≥ 590 nm</b> as reference wavelength.		

**Gentaur Molecular Products  
Voortstraat 49  
1910 Kampenhout, Belgium  
<http://www.gentaur-worldwide.com>**