

RAPID MALARIA PAN-LDH ANTIGEN TEST

A rapid test for the qualitative detection of Malaria pan antigen in human blood sample

Catalog Number:1Q05C2

For in vitro Diagnostic Use Only

INTENDED USE

Rapid Malaria pan-LDH Antigen test is for the rapid qualitative determination of one or more of the known Malaria species; *P. falciparum*, *P. vivax*, *P. ovale*, and/or *P. malariae* by detecting lactate dehydrogenase (LDH) in human blood. This test is an aid in the diagnosis of Malaria infection.

INTRODUCTION

Malaria is a serious parasitic disease characterized by fever, chills, and anemia and is caused by a parasite that is transmitted from one human to another by the bite of infected Anopheles mosquitoes. There are four kinds of malaria that can infect humans: *Plasmodium falciparum*, *Plasmodium vivax*, *Plasmodium ovale*, and *Plasmodium malariae*. In humans, the parasites (called sporozoites) migrate to the liver where they mature and release another form, the merozoites. The disease now occurs in more than 90 countries worldwide, and it is estimated that there are over 500 million clinical cases and 2.7 million malaria-caused deaths per year. At the present, malaria is diagnosed by looking for the parasites in a drop of blood. Blood will be put onto a microscope slide and stained so that the parasites will be visible under a microscope.

PRINCIPLE OF THE TEST

Rapid Malaria pan-LDH Antigen Test contains a membrane strip, which is pre-coated with anti-pan LDH monoclonal antibodies on the test line region of the strip. When a Whole Blood specimen is applied at one end of the membrane and following the application of the assay buffer, it reacts with the colloidal gold-anti-pan LDH antibody that have already been applied to the specimen pad. The mixture then migrates chromatographically towards the other end of the membrane and reacts with the monoclonal antibodies previously placed on the test line region. If the blood contains one or more of the four Malaria species, a colored line will appear in the test line region, showing a positive result. The absence of the colored line in the test region indicates a negative result therefore the whole blood does not contain detectable levels of any of the Malaria species. To serve as a procedural control, a colored line will always appear at the control line region if the test has been performed properly. This control line serves to validate the performance of the test.

MATERIALS PROVIDED

- 1.Rapid Malaria pan Antigen Test Device
- 2.Assay Buffer
- 3.Instructions for Use

MATERIALS NOT PROVIDED

- 1.Calibrated pipette
- 2.Lancet
- 3.Timer

WARNINGS AND PRECAUTIONS

- 1.For professional *in vitro* diagnostic use only. Do not use after expiration date.
- 2.Do not eat, drink or smoke in the area where the specimens or kits are handled.
- 3.Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing and follow the standard procedures for proper disposal of specimens.
- 4.Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being tested.

STORAGE INSTRUCTION

The kit can be stored at room temperature or refrigerated (4-30°C). The test device is stable through the expiration date printed on the sealed pouch. The test device must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

SPECIMEN COLLECTION AND PREPARATION

Collection by venipuncture

- 1) Collect whole blood into a collection tube (containing EDTA, citrate or heparin) by venipuncture.
- 2) If specimens are not immediately tested, they should be refrigerated at 2 ~ 8°C. For storage periods greater than three days, freezing is recommended. They should be brought to room temperature prior to use. Using the specimen after long-term storage of more than three days can cause non-specific reaction.
- 3) When stored at 4 ~ 8°C, the whole blood sample should be used within three days.

Collection using a lancet

- 1) Clean the area to be lanced with an alcohol swab.
- 2) Squeeze the end of the fingertip and pierce with a sterile lancet provided.
- 3) Wipe away the first drop of blood with sterile gauze or cotton.
- 4) Using a 5µL calibrated dropper tube, while gently squeezing the tube, immerse the open end in the blood drop and then gently release the pressure to draw blood into the dropper.

PROCEDURE

Allow test device, buffer, specimen to equilibrate to room temperature(15-30°C) prior to testing.

- 1) Add 5 µl of whole blood into sample well [S1], the small well.
Add two drops (80 µLs) of assay buffer into developer well [S].
- 2) Read the test result in 20 min.

Positive sample may show positive results within a few minutes after the sample is added. However, to confirm a negative result, please wait until 20 minutes.

INTERPRETATION OF RESULTS

1) Positive

The presence of two color bands indicates a positive result for Malaria.


2) Negative

The presence of only one band in the control region of the result window indicates a negative result.

3) Invalid

The test is invalid if the control line does not appear. If this occurs, the test should be repeated using a new device.

(Please refer to the illustrations)



QUALITY CONTROL

1. The control band is an internal reagent and procedural control. It will appear if the test has been performed correctly and the reagents are reactive.
2. Good Laboratory Practice recommends the daily use of control materials to validate the reliability of the device. Control materials which is not provided with this test kit may be commercially available.

LIMITATIONS

- 1) The positive result obtained with Rapid Malaria pan-LDH Antigen Test alone cannot be the final diagnosis of malaria infection. Any positive result must be interpreted in conjunction with the patient clinical history and other laboratory testing results.
- 2) Negative results do not rule out the possibility of malaria exposure or infection.
- 3) As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

ACCURACY

248 patients with suspected uncomplicated malaria were recruited. Blood samples were tested with the rapid malaria pan-LDH antigen test and compared with the gold standard in malaria diagnosis—slide microscopy. The patient sample was assessed as positive or negative by slide microscopy with significant level based on blood parasitemia of 100 parasites/µl of blood determined by expert microscopists/parasitologists.

The test showed true positive in 68 patients, false positive in 7, true negative in 166 and false negative in 7. Sensitivity was 91 % (69/75) and specificity 96 % (166/173). The positive predictive value (PPV) was 91 % (68/75) and the negative predictive value (NPV) 96 % (166/173).

Of all positive samples with rapid malaria pan-LDH test (75), 37(49.3%) was *P.falciparum* alone, 32(42.7%) was mixed infection and 6(8.0%) was positive for non-*P.falciparum* malaria. The lowest level of parasitemia detected by the RDT was 204parasites/µl of blood.

BIBLIOGRAPHY

- 1) Leonard K. Basco, Frederique Marquet, Michael M. Makler, and Jacques Le Bras. : *Plasmodium falciparum* and *Plasmodium vivax*: Lactate Dehydrogenase Activity and its Application for in vitro Drug Susceptibility Assay. Experimental Parasitology 80, 260-271 (1995)
- 2) David L. Vander Jagt, Lucy A. Hunsaker and John E. Heidrich: Partial Purification and Characterization of Lactate Dehydrogenase from *Plasmodium falciparum*. Molecular and Biochemical Parasitology, 4 (1981) 255-264.
- 3) David J. Bzik, Barbara A. Fox and Kenneth Gonyer : Expression of *Plasmodium falciparum* lactate dehydrogenase in *Escherichia coli* Molecular and Biochemical Parasitology, 59(1993) 155-166
- 4) Cameron R. Dunn, Mark J. Banfield, John J. Barker, Christopher W. Highm, Kathleen M. Moreton, Dilek Turgut-Balik, R. Leo Brady and J. John Holbrook. The Structure of lactate dehydrogenase from *Plasmodium falciparum* reveals a new target for anti-malarial design. Nature Structural Biology 3(11)1996, 912-915
- 5) Howard, RJ, et al. The secretion of a Malaria Histidine-rich Protein (Pf HRP II) from Plasmodium falciparum-infected Erythrocytes. J. Cell Biol., (1986)103, 1269-1277
- 6) Rock, EP, et al. Comparative Analysis of the Plasmodium falciparum Histidine-Rich Proteins HRP-I, HRP-II, and HRP-III in Malaria Parasites of Diverse Origin. Parasitol., (1987) 95:209-227
- 7) Parra, ME, et al, Identification of Plasmodium falciparum Histidine-Rich Protein 2 in the plasma of Humans with Malaria. J. Clin. Microbiol., (1991) 29:1629-1634
- 8) Rodriguez-del Valle, M., et al, Detection of Antigens and Antibodies in the Urine of Humans with Plasmodium falciparum Malaria. J. Clin. Microbiol., (1991) 29:1236-1242

