



DATA SHEET

HotStart Taq DNA Polymerase

Store at -20°C

Cat. No.	Description	Concentration	Quantity
G011	HotStart Taq DNA Polymerase	5 U/ μ l	250 U
G039	HotStart Taq DNA Polymerase	5 U/ μ l	1000 U

Product Description

HotStart Taq DNA Polymerase is a chemically modified Taq DNA Polymerase whose enzyme activities can only be activated after 3-5 minutes of incubation at 94°C. This enzyme thus exhibits no polymerase activities before the onset of thermal cycling, preventing non-specific DNA amplification and primer dimer formation. PCR products, amplified up to 6 kb in length with HotStart Taq DNA Polymerase, contain a single base (A) 3' overhang.

Product Components	250 U	1000 U
HotStart Taq DNA Polymerase (5 U/ μ l)	50 μ l	200 μ l
10X PCR buffer, with Mg ²⁺	1 ml	3 ml
25 mM MgSO ₄	1 ml	1 ml

Storage Buffer Components

50 mM Tris-HCl (pH 8.0), 100 mM NaCl, 0.1 mM EDTA, 5 mM DTT, 50 % glycerol and 1.0 % Triton[®]X-100.

Unit Definition

One unit of enzyme catalyzes the incorporation of 10 nmol of deoxyribonucleotides into a polynucleotide fraction in 30 mins at 70°C.

Shipping and Storage

Upon arrival, HotStart Taq DNA Polymerase should be stored at -20°C. Avoid repeated freeze-thaw cycles of all HotStart Taq components to retain maximum performance. All HotStart Taq components are stable for 1 year from the date of shipping if stored and handled properly.

Protocol

The following basic protocol serves as a general guideline and starting point for any PCR amplification. Optimal reaction conditions (incubation times, temperatures, concentration of Taq DNA Polymerase, primers, MgSO₄ and template DNA) may vary and need to be optimized for each specific PCR.

All PCR experiments should be assembled in a nuclease-free environment. In addition, DNA sample preparation, reaction set-up and subsequent reaction(s) should be performed in separate areas to avoid cross contamination.

A negative control reaction (omitting template DNA) should always be performed in tandem with sample PCR to confirm the absence of DNA contamination.

1. Add the following components to a sterile 0.2 ml PCR tube sitting on ice.

Components	Volume	Final Concentration
Template DNA	~100 ng	~2 ng/ μ l
Forward primer (10 μ M)	1 - 2.5 μ l	200 - 500 nM
Reverse primer (10 μ M)	1 - 2.5 μ l	200 - 500 nM
10X PCR buffer, with Mg ²⁺	5 μ l	1X
25 mM MgSO ₄ (optional)*	0 - 3 μ l	1.5 - 3 mM
dNTP Mix (10 mM)	1 μ l	200 μ M
HotStart Taq DNA Polymerase (5 U/ μ l)	0.5 - 1 μ l	2.5 - 5 U
Nuclease-free H ₂ O	up to 50 μ l	-

- * Optimal Mg²⁺ concentration is specific to each DNA template-primer set and can only be determined experimentally.
 - We recommend preparing a mastermix for multiple reactions to minimize reagent loss and enable accurate pipetting.
2. Mix contents of tube and centrifuge briefly.
 3. Incubate tube in a thermal cycler at 94°C for 3 mins to completely denature the template.
 4. Perform 30 - 35 cycles of PCR amplification as follows:
 - Denature:** 94°C for 30 sec
 - Anneal:** 45 - 72°C for 30 sec
 - Extend:** 72°C for 1 min/1 kb template
 5. Incubate for an additional 5 mins at 72°C and maintain the reaction at 4°C. The samples can be stored at -20°C until use.
 6. Analyze the amplification products by agarose gel electrophoresis and visualize by ethidium bromide or SafeView[™] (Cat No. G108) staining. Use appropriate molecular weight standards.

GENTAUR MOLECULAR PRODUCTS
VOORTSTRAAT 49
1910 KAMPENHOUT, BELGIUM