



DATA SHEET

Product name : **Syringe pump can be controlled manually or with a computer. Good for infusion, withdrawal or up to 41 programable pumping phases. Stall detection. Allows pumping of 8 syringes at a time.**

Catalog number : **SP1800**



Similar to SP1000 but holds 2, 4, 6, or 8 syringes (same size) of up to 10 cc.

Infusion rates from 0.001 $\mu\text{L/hr}$ (0.5 μL syringe) to 380 ml/hr (10 cc syringe)

Higher flow rate models available

GENTAUR

Gentaur Molecular Products
Voortstraat 49
1910 Kampenhout, Belgium