

RayBiotech, Inc.

3607 Parkway Lane suite 200 Norcross,GA 30092 Tel: 770-729-2992, 1-888-494-8555

Fax: 770-206-2393

Website: www.raybiotech.com Email: info@raybiotech.com

Certificate of Analysis and Data Sheet

Mouse Anti-Rat CD54 With HRP-conjugated Secondary Antibody

	^^^^^	······	///////
Catalog No.	Species	Isotype	~
SDS-MB-00820	⊗ Human		*
\$20000000000000000000000000000000000000	~~~~~&~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~	~~~~~~~~ ~ ~~~~~~~~~~~~~~~	aaaaaa8

Preparation

Synonyms: ICAM-1

Immunogen: Rat Ax cells (a HEV derived cell line)

Fusion Partners: Spleen cells from immunized BALB/c mice were fused with cells from the PAI

mouse myeloma cell line.

Purification: Purified IgG prepared by affinity chromatography on Protein G from tissue culture

supernatant

Specificity

DS-MB-00820 recognizes the rat CD54 cell surface antigen, also known as intercellular adhesion molecule-1 (ICAM-1), a 90kD adhesion molecule belonging to the immunoglobulin superfamily. CD54 is a cell surface ligand of the lymphocyte integrin, LFA-1 and is known to play an important role in various cell-cell interactions in the immune system. Recent studies suggest that cross-linking of ICAM-1 using clone DS-MB-00820 induces calcium signalling (4). Functionally DS-MB-00820 inhibits homotypic aggregation of PHA blasts.

Formulation

Product Type: Monoclonal Antibody **Product Form:** Purified IgG - liquid

Buffer Solution: Phosphate buffered saline pH7.4 **Preservative Stabilizers:** 0.09% Sodium Azide (NaN₃)

Approx. Protein Concentrations: IgG concentration 1.0 mg/ml

Storage

Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Shelf Life: 18 months from date of dispatch.



RayBiotech, Inc.

3607 Parkway Lane suite 200 Norcross,GA 30092

Tel: 770-729-2992, 1-888-494-8555

Fax: 770-206-2393

Website: www.raybiotech.com Email: info@raybiotech.com

Applications

Table Summary of antibody applications and working conditions

Options Functions	YES	NO	Not determined	Recommended Work dilution or concentration
Flow Cytometry(1)	•			1:25 – 1:50
Immunohistology - Frozen	•			
Immunohistology - Paraffin(2)	•			
Immunohistology - Resin			•	
Immunoprecipitation	•			

Note: Other applications are not tested yet. Optimal dilutions should be determined.

(1) Use 10ul of the suggested working dilution to label 10⁶ cells in 100ul.

(2) This product requires antigen retrieval using heat treatment prior to staining of paraffin sections.

Secondary Antibody Applications

Options Functions	YES	NO	Not determined	Recommended Work dilution or concentration
Immunoassay (ELISA, Western blot)	•			1:5,000-1:10,000

Reference

- 1. Tamatani, T & Miyasaka, M. (1990) Identification of monoclonal antibodies reactive with the rat homolog of ICAM-1, and evidence for a differential involvement of ICAM-1 in the adherence of resting versus activated lymphocytes to high endothelial cells. Int. Immunol. 2: 165-171.
- 2. Tamatani, T. et al. (1991) Molecular mechanisms underlying lymphocyte recirculation. II. Differential regulation of LFA-1 in interaction between lymphocytes and high endothelial cells. Eur. J. Immunol. 21: 855-858.
- 3. Etienne, S. et al. (1998) ICAM-1 signalling pathways associated with Rho activation in microvascular brain endothelial cells. J. Immunol. 161: 5755-5761.
- 4. Etienne, S. et al. (2000) ICAM-1-coupled cytoskeletal rearrangements and transendothelial lymphocyte migration involve intracellular calcium signaling in brain endothelial cell lines. J. Immunol. 165: 3375-3383.



RayBiotech, Inc.

3607 Parkway Lane suite 200 Norcross,GA 30092 Tel: 770-729-2992, 1-888-494-8555

Fax: 770-206-2393

Website: www.raybiotech.com Email: info@raybiotech.com

5. McKechnie, N. M. et al. (2002) Antigenic mimicry: Onchocerca volvulus antigen-specific T cells and ocular inflammation. Invest Ophthalmol Vis Sci. 43:411-8.

6. Azcutia, V. et al. (2010) Inflammation determines the pro-adhesive properties of high extracellular d-glucose in human endothelial cells in vitro and rat microvessels in vivo. PLoS One 5: e10091