

XF-1C Hemoglobin Meter

Usage:

It is used to test anemia.



Feature:

1. Have the function of temperature auto compensation, high stability and measurement accuracy.
2. Advanced photoelectric measurement system greatly reduces the cross pollution.
3. The vertical type has good appearance suitable for regular places such as hospitals, medical labs. The horizontal type is of small size convenient for carrying. There are usually used when medical staff is going on trips for inspection or general examination.
4. The standard calibration memory can automatically save the calibration result for one time. It is very easy and is unnecessary to calibrate repeatedly, which reduces the consumption of calibrating liquid.

GL