

#### RayBiotech, Inc.

3607 Parkway Lane suite 200 Norcross,GA 30092 Tel: 770-729-2992, 1-888-494-8555

Fax: 770-206-2393

Website: www.raybiotech.com Email: info@raybiotech.com

# Certificate of Analysis and Data Sheet

# Recombinant Human B-Cell Leukemia / Lymphoma 2 XL

<	Catalog No.	Source	3
	<sup>3</sup> 228-10105	E.Coli	3
	"~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~		/<

## **Synonyms**

BclXL, Bcl-X(L), Bcl-XL.

#### Introduction

Bcl-XL is a transmembrane protein located in the mitochondrial membranes of cells that are long-lived and postmitotic, such as adult brain cells. It plays a role in the signal transduction pathway of the FAS-Ligand. Bcl-XL is an anti-apoptotic protein which is a member of the Bcl-2 family which are able to form heterodimers, and this is an significant event in the regulation of apoptosis. BCL-XL is involved in the survival of cancer cells.

Bcl-xL is the leading monitor of apoptosis/active cell suicide. Bcl-xL has cell death repressor activity and therefore acts as a survival protein.

# Description

Bcl-XL Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing amino acids 1-210.

The Bcl-XL is purified by proprietary chromatographic techniques.

# Physical Appearance

Sterile Filtered White lyophilized (freeze-dried) powder.

#### **Formulation**

The protein contains 10mM Tris-HCL pH-8, 1mM EDTA and 250mM NaCl.

### Solubility

Suspend BclXL in 100ul of 0.5M Acetic acid, over night at 4°C.

Dilute 10 fold into selected buffer system.Bcl-XL has tendency to form intramolecular disulfide bond, 5mM DTT is recommended in assay buffer. When running SDS-PAGE gel, 10mM DTT is recommended.



### RayBiotech, Inc.

3607 Parkway Lane suite 200 Norcross,GA 30092 Tel: 770-729-2992, 1-888-494-8555

Fax: 770-206-2393

Website: www.raybiotech.com Email: info@raybiotech.com

#### Stability

Lyophilized Bcl-XL although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution Bcl-XL should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

# **Purity**

Greater than 95.0% as determined by SDS-PAGE.