

Recombinant DNA Binding Protein-7 (DBP-7)

CATALOG #:	7603-20	20 µg
	7603-100	100 µg
	7603-1000	1 mg
LOT #:		
SOURCE:	<i>E. coli</i>	
PURITY:	>99% by SDS-PAGE Endotoxin level is <0.1 ng per µg of DBP-7.	
MOL. WEIGHT:	9.44 kDa	
FORM:	Sterile filtered and lyophilized with no additives	

RECONSTITUTION:

Centrifuge the vial prior to opening. Reconstitute in ddH₂O to a concentration of 1.0 mg/ml. Aliquot and store at -20°C for future use. Repeated freeze/thaw cycles should be avoided.

STORAGE CONDITIONS:

The lyophilized DBP-7 is best-stored desiccated below 0°C. Reconstituted protein should be stored in working aliquots at -20°C.


DESCRIPTION:

DNA Binding Protein-7 (DBP-7) is a N-terminal His-tagged recombinant protein and a member of the Sso7d family of small, abundant, non-specific DNA-binding proteins from the hyperthermophilic Archea *Sulfolobus*. The 7-kDa protein from *Sulfolobus* spp. consists of a five stranded, incomplete β-barrel capped at the opening by a C-terminal α-helix; they bind to the minor groove of a DNA duplex via the triple-stranded β-sheet. The topology of the *Sulfolobus* 7-kDa proteins was found to be similar to that of chromatin organization modifier (chromo) domains and eukaryotic SH3 domains, which are involved in protein-protein interactions. *In vitro* studies have shown that DBP-7 promotes the annealing of complementary DNA strands, induces negative supercoiling and chaperones the disassembly and renaturation of protein aggregates in an ATP hydrolysis-dependent manner.

PROTEIN SEQUENCE:

MGSSHHHHHHSSGLVPRGSH**MA**TVKFKYKGEEKEVDISKIKKVVWRVGKMISFTYDEGGGKTGR
GAVSEKDAPKELLQMLEKQKK (Accession# P39476; Gene ID# 1454006).

Note: The bolded methionine above shows the start of DBP-7 sequence (64 aa; 84 aa with His-tag).



Cat. #: 7603-20

Purity of recombinant DNA binding Protein -7 (2 µg) was analyzed by reduced SDS-PAGE electrophoresis,

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptotic Cell Isolation Kit

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- cAMP & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay Kit
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies