

RayBiotech, Inc.

3607 Parkway Lane suite 200 Norcross,GA 30092

Tel: 770-729-2992, 1-888-494-8555

Fax: 770-206-2393

Website: www.raybiotech.com Email: info@raybiotech.com

Certificate of Analysis and Data Sheet

Rabbit Anti Salmonella Species

Catalog No.SpeciesIsotypeMD-05-0418BacteriaRabbit IgG

Preparation

Immunogen: Mixture of S. enteriditis, S. typhimurium, & S. heidelburg

Purification: >95% pure. Protein A chromatography.

Specificity

Polyvalent for Salmonella "O" & "H" antigens. Immunocaptures Salmonellae. Antiserum is not absorbed and does react with related Enterobacteriaceae.

Formulation

Format: Purified, Liquid

Concentration: $4-5 \text{mg/ml} \text{ (OD280nm, E}^{0.1\%} = 1.4)$

Buffer: 0.01M PBS, pH 7.2. This product contains no stabilizing proteins.

Preservative: 0.1% Sodium azide

Storage

Upon receipt, aliquot and store at -20° C. Long term, aliquot and store at -20° C. **Avoid multiple freeze/thaw cycles.**



RayBiotech, Inc.

3607 Parkway Lane suite 200 Norcross,GA 30092

Tel: 770-729-2992, 1-888-494-8555

Fax: 770-206-2393

Website: www.raybiotech.com Email: info@raybiotech.com

Applications

Table Summary of antibody applications and working conditions

Options Functions	YES	NO	Not determined	Recommended Work dilution or concentration
ELISA	•			
Western Blotting			•	
Immunohistology - frozen				
Immunohistology - paraffin				
Immunohistology - resin			•	
Conjugation	•			
Flow Cytometry				
Immunofluorence staining			•	
IFA				

Note: Other applications are not tested yet. Optimal dilutions should be determined by each laboratory for each application.