

#### RayBiotech, Inc.

3607 Parkway Lane suite 200 Norcross,GA 30092

Tel: 770-729-2992, 1-888-494-8555

Fax: 770-206-2393

Website: www.raybiotech.com Email: info@raybiotech.com

## Certificate of Analysis and Data Sheet

# Mouse Anti Salmonella Species (A, B, C, D, E, F, & G Groups)

Catalog No.SpeciesIsotypeMD-05-0419BacteriaIgG2b

#### Preparation

**Source:** Ascites

**Immunogen**: Heat-inactivated, mixture of Salmonella serogroups A, B, C1, C2, and D

**Purification**: >90% pure. Protein A chromatography

### Specificity

Specific for common LPS core of all Salmonella tested. Positive for all Salmonella O-serogroups tested: A, B, C1, C2, D, E1, E3, E4, F, G1, G2. Negative against related bacteria: E.coli (mix), Klebsiella, Citrobacter, Pseudomonas, Yersinia, Shigella, Proteus and Legionella

#### **Formulation**

Format: Purified, Liquid

**Concentration**: Lot specific (OD280nm,  $E^{0.1\%} = 1.4$ )

**Affinity Constant**: Not determined

**Buffer**: 0.01M PBS, pH 7.2. This product contains no stabilizing proteins.

**Preservative**: 0.1% Sodium azide

#### Storage

Upon receipt, aliquot and store at -20°C.

Avoid multiple freeze/thaw cycles.



## RayBiotech, Inc.

3607 Parkway Lane suite 200 Norcross,GA 30092 Tel: 770-729-2992, 1-888-494-8555

Fax: 770-206-2393

Website: www.raybiotech.com Email: info@raybiotech.com

#### **Applications**

Table Summary of antibody applications and working conditions

Options Functions	YES	NO	Not determined	Recommended Work dilution or concentration
ELISA	•			
Western Blotting			•	
Immunohistology - frozen			•	
Immunohistology - paraffin			•	
Immunohistology - resin				
Immunoprecipitation			•	
Flow Cytometry				
Immunofluorence staining			•	
IFA				

Note: Other applications are not tested yet. Optimal dilutions should be determined by each laboratory for each application.