

DATA SHEET

PLX2200 High Frequency remote fluoroscopy X_ray Machine System

Usage

This machine is suitable for human esophagus, chest and stomach, abdomen, limbs and other parts of the work of continuous fluoroscopy diagnosis.

Feature

- 1. The image is clearer for using high quality combined HF HV X-ray generator, high performance image intensifier and high-resolution, high definition CCD.
- 2. With functions of electric up and down movements, electric lateral movement and power rotation, it makes machine move flexibly.
- 3. With the fluoroscopy KV / MA automatic tracking feature, the image brightness and sharpness are at best state.
- 4. Pulse fluoroscope mode reduces patients' absorbing dose and improves the image quality.
- 5. Two control modes-hand brake and foot brake, flexible switching, easy application.
- 6. Programmable LCD with touch screen graphic, simple, elegant.
- 7. Electric iris limited beam, which can be electric rotary, electric opening and closing, electric Oring.
- 8. Microphone makes it convenient for doctors to operate it in other rooms.

