RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

Annexin V Kits & Bulk Reagents Caspase Assay Kits & Reagents

Apoptosis Inducers and Set Apoptosis siRNA Vectors

Nuclear/Cytosol Fractionation Kit Membrane Protein Extraction Kit

Mitochondrial Apoptosis Kits & Reagents Nuclear Apoptosis Kits & Reagents

BioVision

Pimecrolimus (Elidel, SDZ-ASM-981)

CATALOG #: 1593-5 5 mg 1593-25 25 mg

LOT#:

Cell Fractionation System STRUCTURE: Mitochondria/Cytosol Fractionation Kit

> Cytosol/Particulate Rapid Separation Kit Mammalian Cell Extraction Kit

> > FractionPREP Fractionation System Cell Proliferation & Senescence

Quick Cell Proliferation Assay Kit Senescence Detection Kit

High Throughput Apoptosis/Cell Viability Assay Kits

LDH-Cytotoxicity Assay Kit

Bioluminescence Cytotoxicity Assay Kit

Live/Dead Cell Staining Kit

HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits

HAT Colorimetric Assay Kit & Reagents

DNA Damage Quantification Kit

Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

Cell Damage & Repair

cAMP & cGMP Assay Kits

Akt & JNK Activity Assay Kits

Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

Recombinant Adiponectin, Survivin, & Leptin

CETP Activity Assay & Drug Discovery Kits

PLTP Activity Assay & Drug Discovery Kits

Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

siRNA Vectors

Cloning Insert Quick Screening Kit

Mitochondrial & Genomic DNA Isolation Kits

5 Minutes DNA Ligation Kit

20 Minutes Gel Staining/Destaining Kit

β -Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines

Monoclonal and Polyclonal Antibodies

MOLECULAR FORMULA: C₄₃H₆₈CINO₁₁

MOLECULAR WEIGHT: 810.45

CAS 137071-32-0

APPEARANCE: Solid

SOLUBILITY: DMSO (25 mg/ml) or Ethanol (25 mg/ml)

PURITY: >98%

STORAGE: Keep cool and dry, Store at -20°C

DESCRIPTION: Pimecrolimus is a cell-selective inhibitor of inflammatory

> cytokines developed for the treatment of inflammatory skin diseases, such as allergic contact dermatitis, atopic dermatitis, plaque-type psoriasis, and irritant contact dermatitis. It reduces the production of inflammatory cytokines in T cells and mast cells and inhibits the release of preformed

inflammatory mediators from mast cells.

HANDLING: Do not take internally. Wear gloves and mask when handling

the product! Avoid contact by all modes of exposure.

FOR RESEARCH USE ONLY! Not to be used in humans.