BinVision

rev. 01/08

G-CSF Blocking Peptide

50 µg

CATALOG #: 5094RBP-50

AMOUNT:

LOT #:

FORMULATION:

50 μg (0.2 mg/ml) in phosphate buffered saline (PBS), pH 7.2, containing 50% glycerol, 1% BSA and 0.02% thimerosal.

STORAGE CONDITIONS:

Store at -20°C. For long-term storage, aliquot and refreeze at -70°C. Avoid repeated freeze/defrost cycles.

APPLICATION AND USAGE:

The peptide is used for blocking the antibody activity of G-CSF (Cat.#: 5094R-100). It usually blocks the antibody activity completely in Western blot analysis by incubating the peptide with equal volume of antibody for 30-60 minutes at 37°C.

INTERNAL ID#: BV-885-886

FOR RESEARCH USE ONLY! Not to be used on humans.

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents •
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptosis siRNA Vectors

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit ٠

Cell Damage & Repair

.

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits •
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione & Nitric Oxide Fluorometric & Colorimetric Assay Kits Signal Transduction
 - cAMP & cGMP Assay Kits •
 - Akt & JNK Activity Assay Kits
 - Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin •
- CETP Activity Assay & Drug Discovery Kits
- PLTP Activity Assay & Drug Discovery Kits
- **Total Cholesterol Quantification Kit**

Molecular Biology & Reporter Assays

- siRNA Vectors •
- Cloning Insert Quick Screening Kit .
- Mitochondrial & Genomic DNA Isolation Kits .
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit
- β-Galactosidase Staining Kit & Luciferase Reporter Assay Kit

Growth Factors and Cytokines Monoclonal and Polyclonal Antibodies

BioVision

rev. 01/08