

RayBiotech, Inc.

3607 Parkway Lane suite 200 Norcross,GA 30092

Tel: 770-729-2992, 1-888-494-8555

Fax: 770-206-2393

Website: www.raybiotech.com Email: info@raybiotech.com

Certificate of Analysis and Data Sheet

Sheep Anti-BENZOYLECGONINE

With HRP-conjugated Secondary Antibody

Catalog No. DS-PB-00247

Species Chemical

Isotype

Polyclonal IgG

Preparation

Immunogen: Benzoylecgonine conjugated to Keyhole Limpet via the carboxyl group

(Benzoylecgonine-KLH)

Purification: Purified IgG prepared by affinity chromatography on Protein G

Specificity

DS-PB-00247 recognizes benzoylecgonine, a metabolite of cocaine.

It cross-reacts with cocaine but has less than 0.1% reactivity with amphetamine, atropine, codeine, cyclobarbitone, morphine, oxazepam and tetrahydrocannabinol

DS-PB-00247 is suitable for ELISA using Benzoylecgonine-KLH and has been used in FPIA applications

Formulation

Product Form: Purified – liquid

Buffer Solution: Phosphate buffered saline

Preservative & Stabilizers: 0.09 % Sodium Azide

Approx. Protein Concentrations: IgG concentration 5.0 mg/ml

Storage

Store at +4°C or at -20°C if preferred.

This product should be stored undiluted. Storage in frost-free freezers is not recommended. Should this product contain a precipitate we recommend microcentrifugation before use.

Please avoid freeze-thaw cycles.

Shelf Life: 18 months from date of dispatch



RayBiotech, Inc.

3607 Parkway Lane suite 200 Norcross,GA 30092

Tel: 770-729-2992, 1-888-494-8555

Fax: 770-206-2393

Website: www.raybiotech.com Email: info@raybiotech.com

Applications

Table Summary of antibody applications and working conditions

Options Functions	YES	NO	Not determined	Recommended Work dilution or concentration
Immunoassay	•			
Radioimmunoassays	•			

Note: Other applications are not tested yet. Optimal dilutions should be determined by each laboratory for each application.

Secondary Antibody Applications

Options Functions	YES	NO	Not determined	Recommended Work dilution or concentration
Immunoassay (ELISA, Western blot)	•			1:5000-1:10000