



**Innovative**  
R E S E A R C H

**irbiologicals**

**Product:** Cyno Serum

**Lot Number:** Sample

**Volume:** 10ml, 50ml and 100ml Sizes

**Storage:** -20° to -80°

**Description:** Cyno Serum is recovered from whole blood donations from normal healthy Cyno Monkeys of mixed breed. This product is being sold for research and or manufacturing purposes only. It is not to be used in humans or animals. The user assumes all responsibility for its usage and disposal in accordance with all regulations. This material should be handled as if capable of transmitting infectious agents. Please use universal precautions.

**Expiration:** Animal plasma and serum is stable up to three years from the date of receipt when stored at - 20° or colder and not subject to multiple freeze thaw cycles

**Animal Biological References:**

***Determination of the cyclic depsipeptide FK228 in human and mouse plasma by liquid chromatography with mass-spectrometric detection***

X Chen, ER Gardner... - Journal of Chromatography B, 2008 - Elsevier

... Blank mouse plasma was obtained from Innovative Research (Novi, MI). NCI and NCI-Frederick are accredited by AAALAC International and follow the Public Health Service Policy for the Care and Use of Laboratory Animals. ...

***Stability of Caffeic Acid Phenethyl Amide (CAPA) in Rat Plasma***

J Yang, SM Kerwin, PD Bowman... - Biomedical ... - Wiley Online Library

... Heparinized male Sprague–Dawley rat plasma was obtained from Innovative Research (Novi, MI, USA). Water used in the experiments was passed through a Millipore (Billerica, MA, USA) Milli-Q water purification system (18.2 MO). ..

46430 Peary Court Novi, Michigan 48377

ph: 888.660.6866

www.innov-research.com

fx: 248.896.0149

www.irbiologicals.com