

RayBiotech, Inc.

3607 Parkway Lane suite 200 Norcross,GA 30092 Tel: 770-729-2992, 1-888-494-8555

Fax: 770-206-2393

Website: www.raybiotech.com Email: info@raybiotech.com

Certificate of Analysis and Data Sheet

Bovine Heterophile Antigen

Catalog No. Source:

MD-14-0040P Bovine Red Blood cells

Description: Bovine Heterophile Antigen

Format: Purified, Liquid

Purification: Single band by SDS-PAGE (reducing conditions)

Concentration: 2.75mg/ml

Titer: >1:36,000 (ELISA)

Buffer: Deionized water

Preservative: None

Application: Suitable for use in ELISA. Each laboratory should determine an optimum

working titer for use in its particular application. Other applications have not

been tested but use in such assays should not necessarily be excluded.

Inactivation: Not applicable

Storage: Upon receipt, store at -20°C. Avoid multiple freeze/thaw cycles.

References

Clin. Diagn. Virol., (1996, Oct), **7**(1), 17-21.