

#### Trefoil Factor 3 Human E. coli

#### **Product Data Sheet**

Type: Recombinant Cat. No.:

**Source:** E. coli RD172160100 (0.1 mg)

Species: Human

Other names: TFF3, Intestinal trefoil factor, P1.B, ITF,

Polypeptide P1.B, TFI

#### Description

Total 69 AA. MW: 7.82 kDa (calculated). N-terminal His-tag, 10 extra AA (highlighted).

#### Introduction to the Molecule

Trefoil factor 3 (TFF3)belongs together with TFF1 and TFF2 to a small group of mucin-associated peptides. TFF3 contains seven cysteine residues, six of which form disulfide bonds to create a characteristic three-leafed structure. Due to its compact structure, TFF3 is extremely resistant toward acids, proteolytical cleavage or heat degradation. Monomeric form of TFF3 consists of 60 amino acids and has 6.7 kDa, while the dimer (13.1 kDa) consists of 118 amino acids. TFF3 is expressed mainly in gastrointestinal tract, in the mucous cells of the small and large intestine, where it maintains the integrity of mucous layer and in cooperation with mucins protects the gastrointestinal epithelial cells against various injurious agents. However, TFF3 was also detected in salivary glands, posterior pituitary gland and in the inner ear. Secretion of TFF3 is triggered by the presence of certain inflammation mediators and neurotransmitters. Studies showed that oral administration of TFF3 in rats protects gastric mucosa from damage. Over-expression of TFF3 occurs at the sites of damage of the gastrointestinal tract, e.g. peptic ulcer or inflammatory bowel disease. Patients suffering from these diseases have increased levels of TFF3 in serum. TFF3 was reported to be over-expressed also in patients with various neoplasms including intestinal, pancreatic and prostate carcinomas. On the contrary, its expression decreases in thyroid follicular carcinomas. In vitro studies showed that expression of TFF3 is regulated by the level of estrogen in breast cancer cells. Recent study with human and rodent pancreatic islet  $\beta$ -cells has demonstrated that TFF3 over expression increases their proliferation. Both major forms of diabetes involve a decline in islet β-cells mass and their controlled expansion would have great potential utility for treatment of this diseases. Another study with rats has shown that urinary TFF3 protein levels were markedly reduced in response to renal tubular injury, while his levels did not respond to toxicants.

#### Research topic

Energy metabolism and body weight regulation, Oncology, Renal disease, Sepsis

#### **Amino Acid Sequence**

MKHHHHHAS EEYVGLSANQ CAVPAKDRVD CGYPHVTPKE CNNRGCCFDS RIPGVPWCFK PLQEAECTF

## Source

E. coli


## **Purity**

>95%

## SDS-PAGE gel

12% SDS-PAGE separation of Human TFF3

- M.W. marker 14, 21, 31, 45, 66, 97 kDa
- reduced and heated sample, 5µg / lane
- non-reduced and non-heated sample, 5µg / lane



## **Formulation**

Filtered (0.4 µm) and lyophilized from 0.5 mg/ml in 20mM Tris buffer, 150mM NaCl, pH 7.5

#### Reconstitution

Add deionized water to prepare a working stock solution of approximately 0.5 mg/mL and let the lyophilized pellet dissolve completely. Product is not sterile! Please filter the product by an appropriate sterile filter before using it in the cell culture.

### Shipping

At ambient temperature. Upon receipt, store the product at the temperature recommended below.

# Storage, Stability/Shelf Life

Store lyophilized protein at -20°C. Lyophilized protein remains stable until the expiry date when stored at -20°C. Aliquot reconstituted protein to avoid repeated freezing/thawing cycles and store at -80°C for long term storage. Reconstituted protein can be stored at 4°C for a limited period of time; it does not show any change after one week at 4°C.

## **Quality Control Test**

BCA to determine quantity of the protein. SDS PAGE to determine purity of the protein.

# **Applications**

Western blotting

## Note

This product is intended for research use only. The Certificate of Analysis is available on <a href="https://www.biovendor.com">www.biovendor.com</a>					

Gentaur Molecular Products Voortstraat 49 1910 Kampenhout, Belgium http://www.gentaur-worldwide.com